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
IFU Version v4.0.2 (01/06/2025)

INNOVATION SPRINT SRL
Clos Chapelle-aux-Champs 30, bte 1.30.30
1200 Brussels, Belgium

INFORMATIONEN

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Overall
3 studies
422 participants

Tracked Activity
60739 days
128 subjects

Adherence
2804 events from 239 participants
385/624 answered questionnaires

Last week's activity
1 subjects with reported events
1 participants

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Search...

Icon	Name	Organization	Codename	Disease	Patients	Status	Alerts	Actions
	Test Study-I	ORG_1	testirida		0	ONGOING	0	
	Interface	ORG_6	Interface	HIV Remote Patient ...	66	ONGOING	10	
	Symptom Tracker	ORG_1	COVID19		356	ONGOING	285	

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Account Settings

Name: * Therapeutic Area: Codename: *

Protocol: Duration In Months: Start Date: 2021/05/14 End Date:

Primary Language: * Add More Languages: Study Logo: 1

--Select-- Choose File No file chosen

Additional Services

☐ Video-Chat ☐ Multicenter ☐ Interactive Map

MOBILE CONFIGURATION

Study Information Page: Study Logo on App: 1

Leave empty for no Information Page Choose File No file chosen

Pop-up for Study's terms & Privacy 1

Short Term's Description: Terms of Use URL: 1

Leave empty for no Short Terms Description Leave empty for no Terms of Use

Privacy Policy URL: 1

Leave empty for no Privacy Policy

Available Consents

☒ General Terms ☒ Allowing Reporting ☐ Contacted by doctor in case of need

Additional App Widgets

☐ Treatment Reminder ☐ Two-Factor Authentication ☐ Liquid consumption ☐ Sleep Widget

Activity Widget (Available Integrations)

☒ Fitbit Sync ☒ Garmin Sync ☒ Apple Health Sync ☒ Android Sync

SUBJECT DETAIL PAGE CONFIGURATION

☐ First/Last Name ☐ Fiscal Number ☐ Study Dates ☐ Severity

☐ Managed Fitbit ☐ Tracker Description ☐ Mobile Description ☐ Therapeutic Area

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Title	Created On	Type	Status	Actions
Low	23/03/2020	Auto	ACTIVE	[Edit] [Delete]
Medium	23/03/2020	Auto	ACTIVE	[Edit] [Delete]
High	23/03/2020	Auto	ACTIVE	[Edit] [Delete]
Dev	01/04/2020	Custom	ACTIVE	[Edit] [Delete]
Not Synced	29/05/2020	Auto	ACTIVE	[Edit] [Delete]

Abbildung 13 Tags Liste

CREATE DYNAMIC ALERT

Name: *

Status: Active

Portal Alert Message: (You can use {shortid} to include recipient's Short id in your message)

Push Notification Message:

Email List: (Enter emails separated by spaces or commas)

Email Message: (You can use {shortid} to include recipient's Short id in your message)

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ALERTS

Filter by Tags: COVID19 POSITIVE, DOSE 1 COVIDVAX, DOSE 2 COVIDVAX, FEEDBACK, FITBIT TOKEN EXPIRED, HIGH COVID RISK, HIVSRQ SYMPTOM, LOW ADHERENCE, LOW COVID RISK, NOT SYNCED, PAZIENTE FRAGILE, PSYCHIATRIC SUPPORT, PSYCHOLOGICAL SUPPORT

Show: 50

Triggered On (Browser Time)	Subject Id	Alert Message	Tag	Manager
21/04/2021, 14:21:27	20575866	Patient 20575866 needs psychological support	PSYCHOLOGICAL SUPPORT	Giulia Micheli
20/04/2021, 16:30:02	14455901	Subject 14455901 has «Questionario sulla Salute - EQ-5D-3L» as pending		Giulia Micheli
18/04/2021, 16:30:01	04588190	Subject 04588190 has «Questionario sulla Salute - EQ-5D-3L» as pending		Giulia Micheli
18/04/2021, 13:15:19	14455901	Patient 14455901 is identified as FRAIL	PAZIENTE FRAGILE	Giulia Micheli
18/04/2021, 13:12:09	14455901	Subject 14455901 has an HIVSRQ symptom alert	HIVSRQ SYMPTOM	Giulia Micheli
17/04/2021, 19:19:56	03806999	Patient 03806999 needs psychological support	PSYCHOLOGICAL SUPPORT	Giulia Micheli
17/04/2021, 16:30:02	03806999	Subject 03806999 has «Questionario sulla Salute - EQ-5D-3L» as pending		Subject answered the questionnaire.
15/04/2021, 16:30:01	15929649	Subject 15929649 has «Questionario sulla Salute - EQ-5D-3L» as pending		Giulia Micheli
10/04/2021, 09:10:41	11294295	Patient 11294295 needs psychological support	PSYCHOLOGICAL SUPPORT	Cristina Seguti

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Search...

Title	Description	Created On	Status	Actions
test	test	05/03/2021	ACTIVE	[Edit] [Delete]

5 10 All

< 1 >

Abbildung 16 Dynamic List table

Healthentia
Real World Outcomes

Symptom Tracker

DASHBOARD
STUDY
SUBJECTS
USERS
ALERTS
ECONSENT
NOTIFICATIONS
QUESTIONNAIRES
SYSTEM
DEVICES
AUDIT TRAIL

DETAILS DYNAMIC LISTS

Title:
test

Status:
Active

Description:

Patient Criteria

Status:
Enrolled

Gender:
---Select---

Age:

Baseline Date:

Termination Date:

Tags:
Not Synced

Back

Show:
10

Sr. No
Subject ID
Registration Date

No matching records found

Showing 0 to 0 of 0 entries (filtered from 95 total entries)

Figure 17 - Dynamische Listendetails

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Description	Type	Content	Recipients	Sent On (User Time)	Sent On (Your Time)
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...</i>	02/02/2022, 16:24:58
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...</i>	02/02/2022, 14:21:49
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...</i>	31/01/2022, 16:28:38
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...</i>	31/01/2022, 16:27:03
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...</i>	28/01/2022, 20:46:50
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...</i>	28/01/2022, 12:33:32
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...</i>	27/01/2022, 15:43:10
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...</i>	27/01/2022, 14:39:19
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...</i>	21/01/2022, 15:03:14
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...</i>	18/01/2022, 11:54:26

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Description	Type	Content	Status	Actions
Weekly COVID19 Self-Assessment	Weekly	View Content	INACTIVE	Edit Delete
Test	Monthly	View Content	INACTIVE	Edit Delete

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EDIT SCHEDULED NOTIFICATION

Description: Status:

Schedule Type: ☒ Regular ☐ By Baseline Date

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
Title	Codename	Type	No. of Questions	Types	Created On	Status	Actions
Initial Questionnaire	ISPRINT_COVID19	Simple	11	Initial	16/03/2020	ACTIVE	
Body temperature	ISPRINT_FEVER	Simple	1	ADD-EDIT-DELETE	18/03/2020	ACTIVE	
Cough	ISPRINT_COUGH	Simple	2	ADD-EDIT-DELETE	18/03/2020	ACTIVE	
Fatigue	ISPRINT_FATIGUE	Simple	1	ADD-EDIT-DELETE	18/03/2020	ACTIVE	
Difficulty of breath...	ISPRINT_DBREATH	Simple	1	ADD-EDIT-DELETE	26/03/2020	ACTIVE	
Headache	ISPRINT_HEADACHE	Simple	1	ADD-EDIT-DELETE	26/03/2020	ACTIVE	
Muscle Pain	ISPRINT_MUSCLES	Simple	2	ADD-EDIT-DELETE	26/03/2020	ACTIVE	
Abdominal pain	ISPRINT_ABDPAIN	Simple	1	ADD-EDIT-DELETE	26/03/2020	ACTIVE	
Diarrhea	ISPRINT_DIARRHEA	Simple	1	ADD-EDIT-DELETE	26/03/2020	ACTIVE	
Nausea	ISPRINT_NAUSEA	Simple	1	ADD-EDIT-DELETE	26/03/2020	ACTIVE	
Loss of Smell/Taste	ISPRINT_TASTESM	Simple	1	ADD-EDIT-DELETE	31/03/2020	ACTIVE	
Oxygen saturation	ISPRINT_OXYGEN	Simple	1	ADD-EDIT-DELETE	07/04/2020	ACTIVE	
Weekly Questionnaire...	ISPRINT_COVID19_D	Simple	4		08/04/2020	ACTIVE	
Add or take photo	PHOTO	Simple	1	ADD-EDIT-DELETE	24/04/2020	ACTIVE	
Emotional State	EMOSTATUS	Simple	1	ADD-EDIT-DELETE	06/05/2020	INACTIVE	

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The screenshot shows the 'VIEW QUESTIONNAIRE' page in the Healthentia system. The left sidebar contains navigation options: Symptom Tracker, DASHBOARD, STUDY, SUBJECTS, USERS, ALERTS, EDC/ENRICH, NOTIFICATIONS, QUESTIONNAIRES, SYSTEM, DEVICES, and AUDIT TRAIL. The main content area is titled 'VIEW QUESTIONNAIRE' and includes the following fields:

- Title (EN):** Body temperature
- Codename:** ISPRINT_FEVER
- Types:** Add-Event-List
- Status:** Active
- Image:** A placeholder for an image.
- Color:** E55B53
- Description:** A large text area for the questionnaire description.

Below the main form is a 'QUESTION LIST' section. It includes a 'Show' dropdown set to 10, an 'Export' button, and a 'Clear Search Filters' button. The list contains one entry:

No	Code	Question	UI Control	Created On
1	FEVER	What is your body temperature?	Numerical Input	16/03/2020

At the bottom right, there is a 'Cancel' button and a pagination indicator showing 'Showing 1 to 1 of 1 entries'.

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The screenshot shows the 'VIEW QUESTIONNAIRE' page in the Healthentia system for a more complex questionnaire. The left sidebar is the same as in the previous screenshot. The main content area is titled 'VIEW QUESTIONNAIRE' and includes the following fields:

- Title (EN):** Initial Questionnaire
- Codename:** ISPRINT_COVID19
- Types:** Initial
- Status:** Active
- Image:** A placeholder for an image.
- Color:** 808080
- Description:** A large text area for the questionnaire description.

Below the main form is a 'QUESTION LIST' section. It includes a 'Show' dropdown set to 10, an 'Export' button, and a 'Clear Search Filters' button. The list contains five entries:

No	Code	Question	UI Control	Created On
1	BIRTHDT	Enter your Birthdate	Date	13/04/2020
2	GENDER	Define your gender:	Single Choice	16/03/2020
3	DIAGNOSE	Have you been diagnosed with COVID-19?	Single Choice	16/03/2020
4	SYMPTOMS	Do you currently have any of the following symptoms?	Multiple Choice	16/03/2020
5	MEDICIN	Do you usually get any of these medications?	Single Choice	01/04/2020

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ROLES & PERMISSIONS FOR ORG_1

Drag a column header here to group by that column

Search...

Name	Actions
App Manager	Edit Delete
Application	Edit Delete
CRO Assistant	Edit Delete
Investigator	Edit Delete
Organization Admin	Edit Delete
System Administrator	Edit Delete
test	Edit Delete
Viewer	Edit Delete

10 20 50 100 All

CREATE ROLE

Name: *

The Role

Role Access Permissions:

- ☒ All
 - ☒ Dashboard
 - ☒ Export
 - ☒ Subject
 - ☐ Create
 - ☐ Edit
 - ☐ Details
 - ☐ Personal Details
 - ☐ Delete
 - ☐ View Email
 - ☒ Reported Events
 - ☒ Export
 - ☒ Details
 - ☐ Reported Questionnaires
 - ☐ Export
 - ☐ Notifications
 - ☐ Measured Activity Analytics
 - ☐ Consent Actions
 - ☐ Export
 - ☐ Instant Notification
 - ☐ Scheduled Notification
 - ☐ Create
 - ☐ Edit
 - ☐ Delete
 - ☐ Details
 - ☐ ViewContent
 - ☐ Export
 - ☐ Questionnaire
 - ☐ Create
 - ☐ Edit

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The screenshot displays the Healthentia web portal interface. On the left is a purple sidebar with navigation options: Symptom Tracker, DASHBOARD, STUDY, SUBJECTS, USERS (selected), User List, User Invitations, ALERTS, ECONSENT, NOTIFICATIONS, QUESTIONNAIRES, SYSTEM, DEVICES, and AUDIT TRAIL. The main content area is titled 'USER MANAGEMENT' and includes a search bar and a table of users. The table has columns for Email Address, Role, Organization, Created On, and Actions. Below the table are pagination controls showing 5 and 10 items per page, with page 2 of 2 selected. Below the table is the 'CREATE USER' form, which includes fields for First Name, Last Name, Email Address, Password, Confirm Password, Role (dropdown), Time Zone (dropdown), and Sites (text area). At the bottom right of the form are 'Cancel' and 'Submit' buttons.

Email Address	Role	Organization	Created On	Actions
8afe7900-8d69-4030-b546-8c94f2165170@healthentia.com	Application	ORG_1	02/10/2020	[Edit] [More]
bb426815-dbe6-4182-862e-56139ef4d9cd@healthentia.com	System Administrator	ORG_1	01/03/2021	[Edit] [More]
ef9be18f-06ec-4bc1-b551-a8fada21ffe8@healthentia.com	System Administrator	ORG_1	02/10/2020	[Edit] [More]
f0c0b242-19b7-49b3-9b56-cb442ae32b95@healthentia.com	System Administrator	ORG_1	02/10/2020	[Edit] [More]
f6340709-2eb5-48d7-b755-38c5db756c91@healthentia.com	System Administrator	ORG_1	02/10/2020	[Edit] [More]

CREATE USER

First Name: * Last Name: *

Email Address: * Password: Confirm Password:

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Email	Action	Component	Message	Time Stamp
e99be18f-06ec-4bc1-b551-a8fda21...	View	Subject	Subject List on Study: Symptom Trac...	2022/02/03 09:26:32
e99be18f-06ec-4bc1-b551-a8fda21...	Login	User	SUCCESS	2022/02/03 09:26:26
6cabf309-039d-4ade-a5c5-875a763...	View	Medication	Medication List of Subject: Se531aa...	2022/02/03 09:25:39
6cabf309-039d-4ade-a5c5-875a763...	View	PatientQuestionnaire	Subject Questionnaire List of Subjec...	2022/02/03 09:25:39
6cabf309-039d-4ade-a5c5-875a763...	View	EConsent	EConsent List of Subject: Se531aae...	2022/02/03 09:25:39
6cabf309-039d-4ade-a5c5-875a763...	View	Alert	Alert List of Subject: Se531aae-76e4...	2022/02/03 09:25:38
6cabf309-039d-4ade-a5c5-875a763...	View	Conversation	Alert List of Subject: Se531aae-76e4...	2022/02/03 09:25:38
6cabf309-039d-4ade-a5c5-875a763...	View	PatientQuestionnaire	Subject Questionnaire List of Subjec...	2022/02/03 09:25:38
6cabf309-039d-4ade-a5c5-875a763...	View	Subject	Subject Item on Study: Symptom Tra...	2022/02/03 09:25:35
6cabf309-039d-4ade-a5c5-875a763...	View	Subject	Subject List on Study: Symptom Trac...	2022/02/03 09:25:35

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