


INSTRUCTIONS FOR USE

HEALTHCARE PROFESSIONALS

 IFU Version v4.0.2 (01/06/2025)

INNOVATION SPRINT SRL
Clos Chapelle-aux-Champs 30, bte 1.30.30
1200 Brussels, Belgium

ABOUT

The Healthentia Portal is a collection of various modules. Several modules are considered as medical devices and are covered by the CE mark according to the European Regulation 2017/745:



- Subjects-based dashboard: The HEALTHENTIA system offers several subject-level dashboards. There are three broad categories of information to be presented through visualizations related to each subject: subject overview, their reports and their measurements.
- Alerts: The Alerts functionality is using a user interface to set rules and thresholds regarding values collected from various sources like questionnaire answers or IoT devices as an absolute number for a period of time or relative to previous vectors in a period of time. These will then create alerts to clinicians or tag patients.
- Virtual coaching: It is the functionality that allows the user (subject/patient) to interact with an embodied conversational coach for informative and motivational purposes. Users of the HEALTHENTIA Mobile app can interact with the virtual coach through natural language in a controlled manner (I.e., no free text/speech input). The dialogues that are supported by the virtual coach are scripted interactions, which are authored by our clinical experts. Users can navigate through the dialogues by choosing at each step of the conversation between a number of provided reply options. In this way, it is feasible to provide an engaging natural language user interface to the virtual coach, while keeping strict control over which information is provided to the user, or which advice is given.
- Other supportive modules of portal: Study services (questionnaire editor, composite questionnaire, study configurator, multisite study), subject list, study-level dashboard.
- Other supportive modules of backed: scheduler, security & regulatory, management, data handling, API & SDK.

Information

For more information, please visit: <https://healthentia.com>

For any privacy related questions or requests, you can contact: dpo@healthentia.com

For general questions, you can contact: info@healthentia.com

Report of serious incident

Any serious incident that has occurred in relation to the Healthentia Medical Device App, you can contact the manufacturer (Innovation Sprint Srl) in: support@healthentia.com and the authority having jurisdiction in your locale.

Technical support

Healthentia support service is offered at best-effort level.

For any technical support queries, you can contact: support@healthentia.com

Our support team will reply to you within 1 working day.

Paper version of IFU

If you would like a paper version of the instructions for use, please contact the manufacturer at the following e-mail address: info@healthentia.com. It will be delivered within 7 calendar days after receiving the request.

INTENDED PURPOSE

Healthentia is a software intended for: a) the collection and transmission of physiological data including heart rate, blood pressure, oxygen saturation, and weight directly to care providers via automated electronic means in combination with validated IoT devices; b) the visualization (subjects-based dashboards) and the mathematical treatment of data (trends analysis, alerts) related to the monitored chronic disease subject's physiological parameters; c) the transmission of patient's outcomes and outcome scores related to patient's health status, health-affecting factors, health-related quality of life, disease knowledge and adherence to treatment through validated

questionnaires; d) the user (subject/patient) interaction with a conversational virtual coach for informative and motivational purposes, in order to support subject telemonitoring, decision making and virtual coaching.

CLINICAL BENEFITS

The use of Healthentia allows to:

- Provide objective inputs for healthcare professionals to support diagnosis.
- Highlight evolution of physiological parameters by trends analysis of the patient's inputs
- Allow healthcare professionals to provide the same quality of care and safety as the standard of care.

CLINICAL INDICATIONS

Telemonitoring of chronic disease patients (such as heart failure, cancer, COPD, etc.)

CONTRA-INDICATIONS

- Healthentia is not intended for the monitoring of patients in critical situations (operating room, emergency, intensive care).
- Any physical or cognitive condition that, in clinical judgment, would prevent the patient from using Healthentia, e.g. dementia.
- The device should not be used during pregnancy.
- Healthentia has not been tested with paediatric population. Therefore, Healthentia is intended for adult use only (>18 years old) and does not cover use by children.

PATIENTS TARGET GROUP

Chronic disease patients taking part of clinical investigation or a medical treatment

INTENDED USERS

Telemonitored patients and their healthcare professionals

USE ENVIRONMENT AND DURATION

Healthentia can be used as Remote Patient Monitoring solution for patients that are released from hospital and have to follow a certain treatment. Duration depends on the study or intervention or patient's will.

WARNINGS

- The device is not intended to replace the care of a health care professional, including prescription, diagnosis or treatment.
- Consult periodically the Alerts dashboard as it may influence the follow-up of your patients.

PATIENT INFORMATION

Please inform your patient about the following information:

- Healthentia is not monitored in real time by healthcare personnel; it is intended for non-emergency communication only. In case of an emergency, contact your care unit or the emergency number by phone.
- In some special cases, the self-care advice and instructions provided by Healthentia may not be applicable to your situation. If you have questions about your care, symptoms and management, contact the health care team.
- If you have questions about your care, symptoms and management, contact your healthcare team.
- The device is not intended to replace the care of a healthcare professional, including prescription, diagnosis or treatment.
- Telemonitoring does not replace regular check-ups with your doctor. Consult your doctor in case of severe and persistent symptoms.
- In combination with Healthentia, use only measuring devices for which the technical performances (accuracy, precision) have been verified by your healthcare provider as it may affect the outcome of your monitoring.
- Verify the consistency of your clinical data and the good transfer of them to your healthcare provider as the quality of the diagnosis and treatment is partially based on this information.

DISCLAIMER

- The user is responsible of the final assessment of the diagnosis and treatment decided for the patient.
- The device providing information to support diagnostic and therapeutic decision, but the overall clinical context should be taken into consideration before taking any decision.

MINIMUM REQUIREMENTS

Supported hardware

Any personal computer that is using the supported Operating System (OS) and browser – see below- is a supported hardware device.

Operating system/browser requirements

Healthentia Portal requires a stable Internet connection and a compatible device (laptop, desktop, etc). The operating systems supported are:

- Microsoft Windows: version 7 and above
- MacOS: version 10.15 and above
- Linux: LTS distributions released after 2016
- Debian distributions released after 2017
- Fedora distributions released after 2017

This device will need an up-to-date internet browser in order to consult the web application. It is suggested using one of the following browsers that were each tested:

- Google Chrome: version 93 and above
- Mozilla Firefox: version 88 and above
- Safari: version 11 and above
- Microsoft Edge: version 88 and above

Specific issues that come in at Innovation Sprint from the users (post-market surveillance) are analyzed and when a systematic error on one browser is detected, this is added to the backlog and prioritized for the correct version.

INTEGRATION WITH OTHER DEVICES

The accuracy requirements for IoT devices that are connected to Healthentia are listed below.

measurement	Min. accuracy	Medical module
Blood pressure	$\leq 10\text{mmHg}$ (at least 85% probability)	Yes
SpO2	$A_{\text{rms}} \pm 2\text{-}3\%$ of arterial blood gas values	Yes
Heart (RHR, max, ...)	$\pm 10\%$ of the input rate or ± 5 bpm	Yes
Weight	$\pm 0.5\text{-}1.0$ kg	Yes
Physical activity (steps, ...)	n/a	No
Sleep	n/a	No

Healthentia App is compatible with other devices to collect lifestyle information and vital signs. The supported devices intended for use together with Healthentia are:

- Garmin trackers and watches via the Garmin API (Android & iOS)**: Fenix 5 pro, Fenix 6 pro, Forerunner 945, Vivoactive 4, Vivosmart 4, Venu 2S,
- Fitbit trackers and watches via the Fitbit API (Android & iOS)**: Versa, Inspire 2
- iHealth devices via Bluetooth (Android & iOS): Connected Blood Pressure Monitor iHealth Track (KN-550BT)*, Smart body composition scale iHealth Fit (HS2S)**, Smart Pulse Oximeter - iHealth Air (POM3)*
- Polar belt via the Polar API (Android & iOS): H9**

* Device has acceptable accuracy for the intended purpose of Healthentia

** Device does not have acceptable accuracy for the intended purpose of Healthentia and can only be used for measurements that do not require accuracy (e.g. step counter, sleep).

These devices constitute a safe combination and currently there is no device specific information on any known restrictions to combinations. It is the responsibility of the healthcare provider to select the connected devices, among the validated ones listed above, that are appropriate for the intended purpose in terms of measurement accuracy, repeatability and range as it may affect the diagnosis of the patient. Healthentia displays values in international metrics (m, kg, s) regardless of the preference of the patient in his/her IoT device.

SECURITY & PERFORMANCE

Security

Innovation Sprint is committed to protecting the security of data subject's information and takes reasonable precautions to protect it, including protection against unauthorised access, necessary to run the software as intended.

In this context, Innovation Sprint regularly evaluates and reviews technologies, facilities, procedures and potential risks to maintain the security and privacy of users' data, paying particular attention to role-based access and logging of health data access. Moreover, any potential changes to software applications, provided services and/or hardware systems are properly justified, evaluated and registered based on international best practices and standards.

- Avoid using public wi-fi networks to connect to Healthentia when you are about to add and/or manage health data and/or your profile details.
- Do not use jailbroken devices.
- Passwords should be at least eight (8) characters long, contain characters from three of the following four groups: Lowercase letters, Uppercase letters, Numbers (0-9), Special characters.
- Never store your passwords on your devices.
- Change your passwords at least every six months.
- Install an anti-malware program on your devices and update it regularly.
- Make sure that each e-mail you register in Healthentia corresponds to you.
- Ignore and delete messages of doubtful origin and do not follow links that include and refer to websites.
- Do not reveal by phone, e-mail, Internet form, or social media, confidential information such as your username, password.
- Always use the most up-to-date versions of your device operating system and Healthentia App.

Performance

- Provide objective inputs for healthcare professionals to support diagnosis and patient management;
- Highlight evolution of physiological parameters by trends analysis of the patient's inputs;
- Increase patient's adherence to treatment.

TABLE OF CONTENTS

ABOUT.....	2
<i>Intended purpose</i>	<i>2</i>
<i>Clinical benefits</i>	<i>3</i>
<i>Clinical indications.....</i>	<i>3</i>
<i>Contra-indications</i>	<i>3</i>
<i>Patients target group</i>	<i>3</i>
<i>Intended users.....</i>	<i>3</i>
<i>Use environment and duration</i>	<i>3</i>
<i>Warnings.....</i>	<i>3</i>
PATIENT INFORMATION	4
DISCLAIMER.....	4
MINIMUM REQUIREMENTS.....	4
<i>Integration with other devices.....</i>	<i>5</i>
SECURITY & PERFORMANCE	5
<i>Table of Contents</i>	<i>7</i>
1 Getting Started.....	9
1.1 Login/Register as a portal user	9
1.2 Forgot password as a portal user	9
2 STUDY Management	10
2.1 STUDY OVERVIEW.....	10
2.2 Setting up a new STUDY.....	11
3 Dashboard.....	12
3.1 InTERACTIVE MAP.....	12
3.2 BI dashbaord	13
4 Subjects.....	14
4.1 ADD A SUBJECT	14
4.2 SUBJECT RECORD	15

4.3	SUBJECT LIST.....	15
4.4	TAGS & DYNAMIC ALERTS (not a medical module)	16
4.5	ALERTS DASHBOARD	17
4.6	DYNAMIC LISTS.....	17
5	Notifications.....	19
5.1	NOTIFICATION LOG	19
5.2	INSTANT Notifications	19
5.3	Scheduled Notifications	20
6	MESSAGES-TELECONSULTATION.....	21
6.1	MESSAGE INBOX.....	21
6.2	TELECONSULTATION	22
7	ECONSENT.....	23
7.1	eCONSENT	23
8	Questionnaires.....	Error! Bookmark not defined.
8.1	Questionnaire list	25
9	UserS.....	27
9.1	ROLES & PERMISSIONS	27
9.2	User List.....	28
10	AUDIT TRAIL	29
10.1	PORTAL LOG	29
10.2	Application log.....	29

1 GETTING STARTED

HEALTHENTIA can be used both by individuals and by sponsors of clinical trials. Individuals can report and monitor their outcomes (e.g. symptoms) and activity and receive automatic-generated questionnaires for wellbeing. Aggregated data provided after users' consent, can be processed for conducting non-profit research studies. Clinical trial sponsors can operate HEALTHENTIA to collect ePROM/ePREMs under the appropriate regulatory framework (e.g. Good Clinical Practice) to carry out study protocols, often by the use of a PaaS instance of the platform. Individuals, depending on their consent information, may receive invitation to participate in research studies.

This manual is targeting Organization Admins, Study Admin, Investigators and other authorized personnel of the Study Portal.

1.1 LOGIN/REGISTER AS A PORTAL USER

To login user should go to the login page, which is available in <https://saas.healthentia.com> for the SaaS version, or a specified URL for the PaaS version.



Figure 1: Login page




Figure 2 - Registration page

Users receive an invitation from an Organization Admin to register their account and enter in the portal in the allocated study with an allocated role.

1.2 FORGOT PASSWORD AS A PORTAL USER

If you can access to your account, use the forgot password button in the login page.



Copyright © 2022. All Rights Reserved - [Healthentia](#)

Figure 3: Recovery password and email reset

Users receive an email for set up a new password and login to their account.

2 STUDY MANAGEMENT

2.1 STUDY OVERVIEW

After log in you are directed to the Study Overview Page – My Studies, which has a list of your available studies and some overall statistics. You can access your Studies, edit them or create a new one.




Figure 4 - Study Overview

2.2 SETTING UP A NEW STUDY

From this first page, you have the option to create a new study. Depending on the type of study, investigator provides different information. In Figure 5 one can see the available configurations for the study per se but also the mobile App. There is a selection of available widgets to select for the new study and configure the data sources but also several fields to fill in like Name, Code, logo, Languages, therapeutic area, duration and consents.

Healthentia
Real World Outcomes

CREATE STUDY

GENERAL

Account Settings

Name: * Therapeutic Area: Codename: *

Protocol: Duration In Months: Start Date: 2021/06/14 End Date:

Primary Language: * Add More Languages: --Select-- Study Logo: 1

Choose File No file chosen

Additional Services

☐ Video-Chat ☐ Multicenter ☐ Interactive Map

MOBILE CONFIGURATION

Study Information Page:

Leave empty for no Information Page Study Logo on App: 1

Choose File No file chosen

Pop-up for Study's terms & Privacy 1

Short Term's Description:

Leave empty for no Short Terms Description

Terms of Use URL: 1

Leave empty for no Terms of Use

Privacy Policy URL:

Leave empty for no Privacy Policy

Available Consents

☒ General Terms ☒ Allowing Reporting ☐ Contacted by doctor in case of need

Additional App Widgets

☐ Treatment Reminder ☐ Two-Factor Authentication ☐ Liquid consumption ☐ Sleep Widget

Activity Widget (Available integrations)

☒ Fitbit Sync ☒ Garmin Sync ☒ Apple Health Sync ☒ Android Sync

SUBJECT DETAIL PAGE CONFIGURATION

☐ First/Last Name ☐ Fiscal Number ☐ Study Dates ☐ Severity

☐ Managed Fitbit ☐ Tracker Description ☐ Mobile Description ☐ Therapeutic Area

☒ Subject Id ☒ Email ☒ Phone Number ☒ Status

☒ Time Zone ☒ Language ☒ Gender ☒ Weight

☒ Height ☒ Birth Date ☒ Tags ☒ Notes

Cancel Create

Figure 5 - Create New Study

3 DASHBOARD

The user dashboard contains information of the selected study in tabs, e.g. Interactive map, Study Overview, Measured Activities Overview and Outcomes Overview. This is a customizable set of dashboard tabs; they depend on the Outcomes that each trial collects. The map is connected to a Location question, Measured activity requires a connection the mobile sensors (or fitbit, Garmin, etc) and the reported outcomes tab has a self-service mechanisms for questions asked in the questionnaire to be seen in two types of graphs

3.1 INTERACTIVE MAP

At the configuration of the study, the Study Admin can select to enable the interactive map if there are questions regarding location. Portal users can select parameters through several filters and have an overview of subject outcomes at a geographical level, as presented in Figure 6.




Figure 6: Interactive Map

3.2 BI DASHBAORD

Further to the Interactive map, there are other tabs in the Dashboard that can be configured with study Overview statics, like registrations and protocol adherence – Figure 7.




Figure 7: Study registration and compliance overview




Figure 8: Study Measured Activities statistics overview

4 SUBJECTS

4.1 ADD A SUBJECT

Subjects are either imported via an integration with an EDC into Healthentia as inactive subjects that are activated once they have logged in to Healthentia, they can be manually added by the button on the Subject List or invited through the system from the Subject Invitations. The email used for the invitation will be recognized by the system to allow them, at registration on the App to enter to the specific Study. **Not using the same email will not enroll the Subject in the correct Study.**




Figure 9 - Invite Subjects

Similarly, to the Email invitation above, the investigator can invite patients by providing them a code. Subjects can enter this code at the registration, even if they have not received an email. Once they register, they are automatically connected to that specific study only. They can be massively generated to use one per subject and or one can be used for all the subjects for the Study as seen in Figure 10: Create Invitation Codes




Figure 10: Create Invitation Codes

Each Subject has a record of profile details and several other tabs of different information coming from the ePRO or the Activity tracker

4.2 SUBJECT RECORD

In the Subject record, the portal user that has permissions to view or edit the patient details and view ePRO related data and measured activity can have access to the different tabs of information coming from different sources. They can even have the ability to complete questionnaires for them.

Figure 11 - Subject Details

4.3 SUBJECT LIST

The subject's list gathers in a table all patients assigned in the specific study. The list table gives an overview of important elements like dates, Questionnaire adherence, tags, activity status and more.

Email	Subject Id	Registration Date	Qrs Adherence	Tags	Activity Source	Status	Alerts	Actions
008770b1-9470-4993-bb60-...	4N74N	15/04/2020	NaN %	APP TRACKING OPT-OUT	NoTracker	ENROLLED	0	[Edit] [More]
0120b285-3b98-43f5-9b83-...	63PLD	22/10/2020	NaN %	LOW APP TRACKING OPT-OUT	NoTracker	ENROLLED	0	[Edit] [More]
0200764e-c4c2-4777-8d94-...	FR106	25/04/2020	0 %	APP TRACKING OPT-OUT	NoTracker	ENROLLED	1	[Edit] [More]
022c0bd2-897c-40bd-b538-...	Y7B4X	14/04/2020	0 %	HIGH APP TRACKING OPT-OUT	NoTracker	ENROLLED	3	[Edit] [More]
02f1db47-6584-4d33-bae1-...	L4DFD	13/10/2020	NaN %	LOW APP TRACKING OPT-OUT	NoTracker	ENROLLED	0	[Edit] [More]
038dc094-ec78-4d0f-bcb4-6...	XX0I4	30/09/2021	0 %	APP TRACKING OPT-OUT	NoTracker	ENROLLED	1	[Edit] [More]
03042883-bc59-4e1e-a756-...	4KPYB	04/12/2020	NaN %	LOW APP TRACKING OPT-OUT	NoTracker	ENROLLED	0	[Edit] [More]
058d72a6-57fe-49c9-a1d9-...	ZXF3X	15/04/2021	0 %	APP TRACKING OPT-OUT	NoTracker	ENROLLED	1	[Edit] [More]
0630b5a8-b911-43f6-b2f2-2...	FR9E5	30/04/2020	0 %	LOW APP TRACKING OPT-OUT	NoTracker	ENROLLED	2	[Edit] [More]
0733e366-bbdb-4a88-ae2a-...	GVDSJ	08/07/2020	NaN %	LOW APP TRACKING OPT-OUT	NoTracker	ENROLLED	1	[Edit] [More]

Figure 12 Subject list

4.4 TAGS & DYNAMIC ALERTS (not a medical module)

In the section of Tags, we list all the tags available in the study as seen in Figure 13 to be used for grouping the patients in their entry in the system or depending on their answers in the questionnaires. These tags can be then attached to an Alert seen in Figure 14 - Dynamic Alerts List. All Alerts are featured in a table as seen in Figure 14. As in all sections of the portal all data are exportable to an excel file.

Title	Created On	Type	Status	Actions
Low	23/03/2020	Auto	ACTIVE	Edit Delete
Medium	23/03/2020	Auto	ACTIVE	Edit Delete
High	23/03/2020	Auto	ACTIVE	Edit Delete
Dev	01/04/2020	Custom	ACTIVE	Edit Delete
Not Synced	29/05/2020	Auto	ACTIVE	Edit Delete

Figure 13 Tags List

CREATE DYNAMIC ALERT

Name:

Status:

Portal Alert Message: (You can use {shortid} to include recipient's Short id in your message)

Push Notification Message:

Email List: (Enter emails separated by spaces or commas)

Email Message: (You can use {shortid} to include recipient's Short id in your message)

API Endpoint to Call:

API Alert Codename:

Tags to Apply:

Figure 14 - Dynamic Alerts List

4.5 ALERTS DASHBOARD

The Alerts are shown in a separate Dashboard where the portal user has an overview of what tags are allocated and whether a patients requires an immediate action or has low adherence in his reporting of questionnaires.

ALERTS

Filter by Tags: COVID19 POSITIVE, DOSE 1 COVIDVAX, DOSE 2 COVIDVAX, FEEDBACK, FITBIT TOKEN EXPIRED, HIGH COVID RISK, HIVSRQ SYMPTOM, LOW ADHERENCE, LOW COVID RISK, NOT SYNCED, PAZIENTE FRAGILE, PSYCHIATRIC REDDIT, PSYCHOLOGICAL SUPPORT

Show: 50

Triggered On (Browser Time)	Subject Id	Alert Message	Tag	Manager
21/04/2021, 14:21:27	20575866	Patient 20575866 needs psychological support	PSYCHOLOGICAL SUPPORT	Giulia Micheli
20/04/2021, 16:30:02	14455901	Subject 14455901 has «Questionario sulla Salute - EQ-5D-3L» as pending		Giulia Micheli
18/04/2021, 16:30:01	04588190	Subject 04588190 has «Questionario sulla Salute - EQ-5D-3L» as pending		Giulia Micheli
18/04/2021, 13:15:19	14455901	Patient 14455901 is identified as FRAIL	PAZIENTE FRAGILE	Giulia Micheli
18/04/2021, 13:12:09	14455901	Subject 14455901 has an HIVSRQ symptom alert	HIVSRQ SYMPTOM	Giulia Micheli
17/04/2021, 19:19:56	03806999	Patient 03806999 needs psychological support	PSYCHOLOGICAL SUPPORT	Giulia Micheli
17/04/2021, 16:30:02	03806999	Subject 03806999 has «Questionario sulla Salute - EQ-5D-3L» as pending		Subject answered the questionnaire.
15/04/2021, 16:30:01	15929649	Subject 15929649 has «Questionario sulla Salute - EQ-5D-3L» as pending		Giulia Micheli
10/04/2021, 09:10:41	11294295	Patient 11294295 needs psychological support	PSYCHOLOGICAL SUPPORT	Cristina Seguti

Figure 15 - Alerts Dashboard

4.6 DYNAMIC LISTS

The Dynamic lists are used to group patients in a dynamic way. We create grouping requirements and patients are added automatically in a list when they meet the criteria. The attributes of patients that can be used as criteria for filtering the patient list are status, screening, termination or treatment switch date and even using tags that are explained below. As in all sections of the portal all data are exportable to an excel file, as presented in Figure 16.

DYNAMIC LISTS

Drag a column header here to group by that column

Search...

Title	Description	Created On	Status	Actions
test	test	05/03/2021	ACTIVE	[Export] [Delete]

5 10 All

< 1 >

Figure 16 Dynamic List table

Healthentia
Real World Outcomes

Symptom Tracker

DASHBOARD
STUDY
SUBJECTS
USERS
ALERTS
ECONSENT
NOTIFICATIONS
QUESTIONNAIRES
SYSTEM
DEVICES
AUDIT TRAIL

DETAILS DYNAMIC LISTS

Title:
test
Status:
Active

Description:

Patient Criteria

Status:
Enrolled
Gender:
---Select---
Age:

Baseline Date:
Termination Date:
Tags:
Not Synced

Show:
10

Sr. No	Subject ID	Registration Date
No matching records found		

Showing 0 to 0 of 0 entries (filtered from 95 total entries)

Figure 17 - Dynamic list details

5 NOTIFICATIONS

5.1 NOTIFICATION LOG

In the section of Notifications, we can see all communication sent by the system to patients, we can create instant notifications or scheduled. All notifications are listed by date in the Notification List and you have an overlook of the content of these notifications and if their type. As in all sections of the portal all data are exportable to an excel file, as presented in Figure 18.

Description	Type	Content	Recipients	Sent On (User Time)	Sent On (Your Time)
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...	02/02/2022, 16:24:58
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...	02/02/2022, 14:21:49
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...	31/01/2022, 16:28:38
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...	31/01/2022, 16:27:03
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...	28/01/2022, 20:46:50
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...	28/01/2022, 12:33:32
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...	27/01/2022, 15:43:10
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...	27/01/2022, 14:39:19
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...	21/01/2022, 15:03:14
<i>Instant Notification</i>	Instant	View Content	View Recipients (1)	<i>Based on each user's time...	18/01/2022, 11:54:26

Figure 18 - Notification Log

5.2 INSTANT NOTIFICATIONS

Instant Notifications are the ad hoc messages that the investigator can send to a patient that hold a simple message or accompany a questionnaire. You can select individual recipients for the notification from the full subject list or filter by using tags or Dynamic Lists as seen in Figure 19.

Figure 19 - Instant Notification

5.3 SCHEDULED NOTIFICATIONS

Scheduled notifications are the ones that are created in the beginning of the study and are used to schedule the questionnaires to be sent to a group of recipients based on a filtering by tags or dynamic lists or to the whole list.

The scheduling of a notification can be done based on a regular Date with a Daily, weekly or Monthly interval sent once at a specific time or recurring. The system also allows the coordinator to send notifications based on the protocol's specific dates and X days before or after these days. You can see the Notifications List at Figure 20 and the New Scheduling Notification at Figure 21 - New Scheduled Notification.

Description	Type	Content	Status	Actions
Weekly COVID19 Self-Assessment	Weekly	View Content	INACTIVE	Edit Delete
Test	Monthly	View Content	INACTIVE	Edit Delete

Figure 20 - Scheduled Notifications

EDIT SCHEDULED NOTIFICATION

Description: Status:

Schedule Type: ☒ Regular ☐ By Baseline Date

Frequency: Day of Month: End repeat:

Recipient Tags:

Questionnaires:

Message (Individual Notification):

Figure 21 - New Scheduled Notification

6 MESSAGES-TELECONSULTATION

6.1 MESSAGE INBOX

End-users of the mobile app if the feature is enabled they can send a Message to the Investigator/study assistant and start a conversation. If needed the investigator can start a teleconsultation video link that will be valid for 15 minutes. See inbox in Figure 22 and the conversation in Figure 23.

Figure 23.




Figure 22 - Message Inbox




Figure 23 - Message conversation

6.2 TELECONSULTATION

If needed, the investigator can start a teleconsultation video link that will be valid for 15 minutes, as shown in the Figure below.




Figure 24 – Teleconsultation

7 ECONSENT

7.1 ECONSENT

If the trial has an active eConsent, the eConsent is triggered by the Healthentia app after the user has registered and agreed with the Terms of the Healthentia application. Each eConsent form must have a unique Version and a unique Title in the context of this study or site (in case of multicenter trials). When creating or editing an eConsent form, users with the appropriate permission shall be able to: (1) add a plain text field to write the "Terms" of the eConsent Form; (2) add a plain text field to write the "Introduction" of the eConsent form; (3) upload a url from a video streaming provider to the eConsent form; and (4) define a Status. An eConsent setup has three status: "draft", "active" and "inactive". While in "draft" mode, the eConsent is editable. Once the investigator changes its status to "active", the eConsent is locked and no longer editable. The "inactive" status is given to eConsent forms that are no longer active but that have been signed at least by one study participant. See Figure 25 and Figure 26.



Figure 25 - eConsent list




Figure 26 - eConsent configuration

The process of signing an eConsent form can only be completed after a mobile app user has verified his identity providing a code received via email or mobile phone. Once the participant has signed the eConsent, a pdf document is created with the name, date and signature of the participant on top of the Terms plain text. The pdf document created after the participant signed the eConsent is saved in the Healthentia Patient documents with type "consent". Once the participant has signed the eConsent, the tag of the eConsent in the subject list changes from "not started" to "requires signature".

Healthentia
Real World Outcomes

Symptom Tracker

DASHBOARD

STUDY

SUBJECTS

USERS

ALERTS

ECONSENT

eConsent Setup

Patient Status

NOTIFICATIONS

QUESTIONNAIRES

SYSTEM

DEVICES

AUDIT TRAIL

SUBJECT ECONSENT LIST

Drag a column header here to group by that column

Search...

10

20

Sr. No	Email	Short Id	Consent Info	Version	Consent Date	Comprehension	Status	File
1	008770b1-9470-499...	4N74N	ECONSENT	1			NOT REQUIRED	
2	0120b285-3b98-43f...	63PLD	ECONSENT	1			NOT REQUIRED	
3	0200764e-c4c2-477...	FR106	ECONSENT	1			NOT REQUIRED	
4	022c0bd2-897c-40b...	Y7B4X	ECONSENT	1			NOT REQUIRED	
5	02f1db47-6584-4d3...	L4DFD	ECONSENT	1			NOT REQUIRED	
6	0386c094-ec78-460...	XXOI4	ECONSENT	1			NOT REQUIRED	
7	03be2883-bc59-4e1...	4KPYB	ECONSENT	1			NOT REQUIRED	
8	058672a6-57fe-49c...	ZXF3X	ECONSENT	1			NOT REQUIRED	
9	0636b5a8-b911-43f...	FR9E5	ECONSENT	1			NOT REQUIRED	
10	0733e366-bbdb-4a...	GVD5J	ECONSENT	1			NOT REQUIRED	

1

2

3

4

5

...

36

Figure 27 – Subject eConsent list

8 QUESTIONNAIRES

8.1 QUESTIONNAIRE LIST

The portal allows the study coordinators to create a questionnaire with different question types but also allows the questionnaire to be uploaded from a specified excel file. These questionnaires are then sent to the patient through a notification. In the Questionnaire management one can see the list of available questionnaires in the study and the number of questions that each has as well as their status and date created.

A list of qualified questionnaires is made available. A qualified questionnaire, specifically a Patient-Reported Outcome Measure (PROM), in the context of a clinical research system, refers to a standardized, validated instrument used to assess health outcomes directly from patients. Such qualified questionnaires include standardized set of questions and be validated through rigorous scientific processes to ensure reliability, validity, and responsiveness to change. Qualified questionnaires are identified by a tag in the Questionnaire Management tool. This subset of questionnaires is part of the medical modules of Healthentia and are subject to the CE mark. This is also described on the website and a list is available in the portal of the healthcare professionals.

As in all sections of the portal all questionnaires can be exportable to an excel file.

Title	CodeName	Type	No. of Questions	Types	Created On	Status	Actions
Initial Questionnaire	ISPRINT_COVID19	Simple	11	INITIAL	16/03/2020	ACTIVE	
Body temperature	ISPRINT_FEVER	Simple	1	ADD-EVENT-LIST	18/03/2020	ACTIVE	
Cough	ISPRINT_COUGH	Simple	2	ADD-EVENT-LIST	18/03/2020	ACTIVE	
Fatigue	ISPRINT_FATIGUE	Simple	1	ADD-EVENT-LIST	18/03/2020	ACTIVE	
Difficulty of breath...	ISPRINT_DIBREATH	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	
Headache	ISPRINT_HEADACHE	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	
Muscle Pain	ISPRINT_MUSCLES	Simple	2	ADD-EVENT-LIST	26/03/2020	ACTIVE	
Abdominal pain	ISPRINT_ABDRAIN	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	
Diarrhea	ISPRINT_DIARRHEA	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	
Nausea	ISPRINT_NAUSEA	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	
Loss of Smell/Taste	ISPRINT_TASTESM	Simple	1	ADD-EVENT-LIST	31/03/2020	ACTIVE	
Oxygen saturation	ISPRINT_OXYGEN	Simple	1	ADD-EVENT-LIST	07/04/2020	ACTIVE	
Weekly Questionnaire...	ISPRINT_COVID19_D	Simple	4	ADD-EVENT-LIST	08/04/2020	ACTIVE	
Add or take photo	PHOTO	Simple	1	ADD-EVENT-LIST	24/04/2020	ACTIVE	
Emotional State	EMOSTATUS	Simple	1	ADD-EVENT-LIST	06/05/2020	INACTIVE	

Figure 28 Questionnaire Management

When creating a new one you define the title and questionnaire type (select from a list created by the admin) and upload a characteristic image for it if you want it to be shown on the app and start adding the questions one by one. In the question options you can choose from a variety of UI controls like if the question is a single, multiple, entry text or slider bar as presented in Figure 29.

The screenshot shows the 'VIEW QUESTIONNAIRE' page in the Healthentia web portal. The left sidebar contains navigation options: Symptom Tracker, DASHBOARD, STUDY, SUBJECTS, USERS, ALERTS, ECOSYSTEM, NOTIFICATIONS, QUESTIONNAIRES, SYSTEM, DEVICES, and AUDIT TRAIL. The main content area is titled 'VIEW QUESTIONNAIRE' and includes the following fields:

- Title (EN):** Body temperature
- Codename:** ISPRINT_FEVER
- Types:** Add-Event-List
- Status:** Active
- Image:** A placeholder image with a red border.
- Color:** A red color swatch with the code E55B53.
- Description:** A large text area for the questionnaire description.

Below these fields is a 'QUESTION LIST' section. It includes a 'Show: 10' dropdown, an 'Export' button, and a 'Clear Search Filters' button. The list table has columns: No, Code, Question, UI Control, and Created On. The table contains one entry:

No	Code	Question	UI Control	Created On
1	FEVER	What is your body temperature?	Numerical Input	16/03/2020

At the bottom right, there is a 'Cancel' button and a pagination indicator showing '1' of 1 entries.

Figure 29 - Simple Questionnaire Details

The screenshot shows the 'VIEW QUESTIONNAIRE' page in the Healthentia web portal for a composite questionnaire. The left sidebar is the same as in Figure 29. The main content area is titled 'VIEW QUESTIONNAIRE' and includes the following fields:

- Title (EN):** Initial Questionnaire
- Codename:** ISPRINT_COVID19
- Types:** Initial
- Status:** Active
- Image:** A placeholder image with a red border.
- Color:** A dark gray color swatch with the code 808080.
- Description:** A large text area for the questionnaire description.

Below these fields is a 'QUESTION LIST' section. It includes a 'Show: 10' dropdown, an 'Export' button, and a 'Clear Search Filters' button. The list table has columns: No, Code, Question, UI Control, and Created On. The table contains five entries:

No	Code	Question	UI Control	Created On
1	BIRTHDT	Enter your Birthdate	Date	13/04/2020
2	GENDER	Define your gender:	Single Choice	16/03/2020
3	DIAGNOSE	Have you been diagnosed with COVID-19?	Single Choice	16/03/2020
4	SYMPTOMS	Do you currently have any of the following symptoms?	Multiple Choice	16/03/2020
5	MEDICIN	Do you usually get any of these medications?	Single Choice	01/04/2020

Figure 30 - Composite Questionnaire

9 USERS

9.1 ROLES & PERMISSIONS

This section is used by the admin to create the main roles and their permissions in the portal sections. Permissions can be configured for viewing, editing and delete, as depicted in Figure 31.

The figure consists of two screenshots of the Healthentia web portal interface, showing the 'ROLES & PERMISSIONS' section.

Top Screenshot: ROLES & PERMISSIONS FOR ORG_1

The interface shows a sidebar with navigation options: Symptom Tracker, DASHBOARD, STUDY, SUBJECTS, USERS, ALERTS, ECONSENT, NOTIFICATIONS, QUESTIONNAIRES, SYSTEM, DEVICES, and AUDIT TRAIL. The main content area is titled 'ROLES & PERMISSIONS FOR ORG_1'. It includes a search bar and a table of roles.

Name	Actions
App Manager	View Edit Delete
Application	View Edit Delete
CRO Assistant	View Edit Delete
Investigator	View Edit Delete
Organization Admin	View Edit Delete
System Administrator	View Edit Delete
test	View Edit Delete
Viewer	View Edit Delete

At the bottom of the table, there are pagination options: 10, 20, 50, 100, and All.

Bottom Screenshot: CREATE ROLE


The interface shows the 'CREATE ROLE' form. It includes a 'Name' field with the placeholder text 'The Role'. Below the name field, there is a section for 'Role Access Permissions' with a tree structure of permissions.

- ☒ All
 - ☒ Dashboard
 - ☒ Export
 - ☒ Subject
 - ☐ Create
 - ☐ Edit
 - ☐ Details
 - ☐ Personal Details
 - ☐ Delete
 - ☐ View Email
 - ☒ Reported Events
 - ☒ Export
 - ☒ Details
 - ☐ Reported Questionnaires
 - ☐ Export
 - ☐ Notifications
 - ☐ Measured Activity Analytics
 - ☐ Consent Actions
 - ☐ Export
 - ☐ Export
 - ☐ Instant Notification
 - ☐ Scheduled Notification
 - ☐ Create
 - ☐ Edit
 - ☐ Delete
 - ☐ Details
 - ☐ ViewContent
 - ☐ Export
 - ☐ Questionnaire
 - ☐ Create
 - ☐ Edit

Figure 31 Role Management and Add new Role

9.2 USER LIST

In the users list you can find all the portal users that are assigned in the study. The study admin can create or invite new users or assign existing ones to a specific study with a permission role. The user provides the email of an investigator that triggers an auto generated email to that email. Upon receipt, investigators need to confirm the email and follow a registration. This list can be then exported to an excel, as depicted in Figure 32.



USER MANAGEMENT

Drag a column header here to group by that column

Search...

Email Address	Role	Organization	Created On	Actions
6afe7900-8d69-4030-b546-8c94f2165170@healthentia.com	Application	ORG_1	02/10/2020	[Edit] [More]
bb426815-dbe6-4182-862e-56139ef4d9cd@healthentia.com	System Administrator	ORG_1	01/03/2021	[Edit] [More]
ef9be18f-06ec-4bc1-b551-a8fada21ffe8@healthentia.com	System Administrator	ORG_1	02/10/2020	[Edit] [More]
f0c0b242-19b7-49b3-9b56-cb442ae32b95@healthentia.com	System Administrator	ORG_1	02/10/2020	[Edit] [More]
f8340709-2eb5-48d7-b755-38c5db756c91@healthentia.com	System Administrator	ORG_1	02/10/2020	[Edit] [More]

5 10 < 1 2 >

CREATE USER

First Name: * Last Name: *

Email Address: * Password: * Confirm Password: *

Role: * Time Zone: *

Sites:

Cancel Submit

Figure 32 User Management and User Creation

10 AUDIT TRAIL

10.1 PORTAL LOG

In the section of Log Management all actions that users do from viewing to editing or deleting. The log list of actions can be filtered by date from the top of the page, while each column on the table allows the search for a specific result. You can also sort by alphabetic or numeric order per column. When viewing the log details you can get more information on the reported action. As in all sections of the portal all data are exportable to an excel file.

PORTAL LOG

Drag a column header here to group by that column

Email	Action	Component	Message	Time Stamp
ef9be18f-06ec-4bc1-b551-a8fada21...	View	Subject	Subject List on Study: Symptom Trac...	2022/02/03 09:26:32
ef9be18f-06ec-4bc1-b551-a8fada21...	Login	User	SUCCESS	2022/02/03 09:26:26
6cabf309-039d-4ade-a5c5-875a763...	View	Medication	Medication List of Subject: Se531aa...	2022/02/03 09:25:39
6cabf309-039d-4ade-a5c5-875a763...	View	PatientQuestionnaire	Subject Questionnaire List of Subjec...	2022/02/03 09:25:39
6cabf309-039d-4ade-a5c5-875a763...	View	EConsent	EConsent List of Subject: Se531aae...	2022/02/03 09:25:39
6cabf309-039d-4ade-a5c5-875a763...	View	Alert	Alert List of Subject: Se531aae-76e4...	2022/02/03 09:25:38
6cabf309-039d-4ade-a5c5-875a763...	View	Conversation	Alert List of Subject: Se531aae-76e4...	2022/02/03 09:25:38
6cabf309-039d-4ade-a5c5-875a763...	View	PatientQuestionnaire	Subject Questionnaire List of Subjec...	2022/02/03 09:25:38
6cabf309-039d-4ade-a5c5-875a763...	View	Subject	Subject Item on Study: Symptom Tra...	2022/02/03 09:25:35
6cabf309-039d-4ade-a5c5-875a763...	View	Subject	Subject List on Study: Symptom Trac...	2022/02/03 09:25:35

10 20

Figure 33 Log Management & Details

10.2 APPLICATION LOG

The application log, provides information about the actions of the subjects, e.g. their participation in questionnaires, as shown in Figure

APPLICATION LOG

Drag a column header here to group by that column

Subject ID	Description	Log Type	Component	Time Stamp
O31ES	FAILED	FailedLogin	Subject	2022/02/03 09:27:34
O31ES	LOG_OUT	Logout	Subject	2022/02/03 09:24:00
O31ES	SUCCESS (NOT_CONFIRMED_EMAIL) Login	Login	Subject	2022/02/03 08:26:38
JAQYU	LOG_OUT	Logout	Subject	2022/02/03 08:26:22
JAQYU	SUCCESS (NOT_CONFIRMED_EMAIL) Login	Login	Subject	2022/02/03 08:25:28
O31ES	LOG_OUT	Logout	Subject	2022/02/03 08:25:17
DMT6E	SUCCESS (NOT_CONFIRMED_EMAIL) Login	Login	Subject	2022/02/02 03:01:34
DMT6E	SUCCESS (NOT_CONFIRMED_EMAIL) Login	Login	Subject	2022/02/02 02:54:14
O31ES	SUCCESS (NOT_CONFIRMED_EMAIL) Login	Login	Subject	2022/02/02 02:33:03
DMT6E	SUCCESS (NOT_CONFIRMED_EMAIL) Login	Login	Subject	2022/02/02 01:25:24

10 20

Figure 34 Application Log Detail