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Copyright © 2022. All Rights Reserved - Healthentia

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Icon	Name	Organization	Codename	Disease	Patients	Status	Alerts	Actions
	Test Study-1	ORG_1	testirida		0	ONGOING	0	
	Interface	ORG_6	Interface	HIV Remote Patient ...	66	ONGOING	10	
	Symptom Tracker	ORG_1	COVID19		356	ONGOING	285	

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**Account Settings**

Name: \* Therapeutic Area: Codename: \*


Protocol: Duration In Months: Start Date: 2021/05/14 End Date:

Primary Language: \* Add More Languages: Study Logo:  Choose File No file chosen

**Additional Services**

Video-Chat  Multicenter  Interactive Map

**MOBILE CONFIGURATION**

Study Information Page: Leave empty for no information Page Study Logo on App:  Choose File No file chosen

Pop-up for Study's terms & Privacy  Short Term's Description: Leave empty for no Short Terms Description Terms of Use URL:  Leave empty for no Terms of Use Privacy Policy URL: Leave empty for no Privacy Policy

**Available Consents**

General Terms  Allowing Reporting  Contacted by doctor in case of need

**Additional App Widgets**

Treatment Reminder  Two-Factor Authentication  Liquid consumption  Sleep Widget

**Activity Widget (Available Integrations)**

Fitbit Sync  Garmin Sync  Apple Health Sync  Android Sync

**SUBJECT DETAIL PAGE CONFIGURATION**

First/Last Name  Fiscal Number  Study Dates  Severity

Managed Fitbit  Tracker Description  Mobile Description  Therapeutic Area

Subject Id  Email  Phone Number  Status

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

The screenshot shows the 'TAG MANAGEMENT' section of the Healthentia web portal. The left sidebar contains navigation options: Symptom Tracker, DASHBOARD, STUDY, SUBJECTS (with sub-options: Subject List, Tag List, Dynamic Lists, Custom Fields, Subject Invitations, Codes), and Codes. The main content area displays a table of tags with the following columns: Title, Created On, Type, Status, and Actions. The table contains five rows of data:

Title	Created On	Type	Status	Actions
Low	23/03/2020	Auto	ACTIVE	[Edit] [Delete]
Medium	23/03/2020	Auto	ACTIVE	[Edit] [Delete]
High	23/03/2020	Auto	ACTIVE	[Edit] [Delete]
Dev	01/04/2020	Custom	ACTIVE	[Edit] [Delete]
Not Synced	29/05/2020	Auto	ACTIVE	[Edit] [Delete]

At the bottom of the table, there are pagination controls showing '5' items per page, '10' total items, and 'All' filter options. The interface also includes a search bar and various utility icons at the top right.

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The screenshot shows the 'CREATE DYNAMIC ALERT' form in the Healthentia web portal. The left sidebar is identical to the previous screenshot. The main content area contains the following fields:

- Name:** A text input field.
- Status:** A dropdown menu with 'Active' selected.
- Portal Alert Message:** A text input field with a note: '(You can use [shortid] to include recipient's Short id in your message)'. Below it is a small text input for a placeholder.
- Push Notification Message:** A text input field.
- Email List:** A text input field with a note: '(Enter emails separated by spaces or commas)'. Below it is a small text input for a placeholder.
- Email Message:** A text input field with a note: '(You can use [shortid] to include recipient's Short id in your message)'. Below it is a small text input for a placeholder.
- API Endpoint to Call:** A text input field.
- API Alert Codename:** A text input field.
- Tags to Apply:** A text input field.

At the bottom right, there are 'Cancel' and 'Create' buttons.

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**ALERTS**

Filter by Tags: COVID19 POSITIVE, DOSE 1 COVIDVAX, DOSE 2 COVIDVAX, FEEDBACK, FITBIT TOKEN EXPIRED, HIGH COVID RISK, HIVSRQ SYMPTOM, LOW ADHERENCE, LOW COVID RISK, NOT SYNCED, PAZIENTE FRAGILE, PSYCHIATRIC SUPPORT, PSYCHOLOGICAL SUPPORT

Show: 50

Triggered On (Browser Time)	Subject Id	Alert Message	Tag	Manager
21/04/2021, 14:21:27	20575866	Patient 20575866 needs psychological support	PSYCHOLOGICAL SUPPORT	Giulia Micheli
20/04/2021, 16:30:02	14455901	Subject 14455901 has «Questionario sulla Salute - EQ-5D-3L» as pending		Giulia Micheli
18/04/2021, 16:30:01	04588190	Subject 04588190 has «Questionario sulla Salute - EQ-5D-3L» as pending		Giulia Micheli
18/04/2021, 13:15:19	14455901	Patient 14455901 is identified as FRAIL	PAZIENTE FRAGILE	Giulia Micheli
18/04/2021, 13:12:09	14455901	Subject 14455901 has an HIVSRQ symptom alert	HIVSRQ SYMPTOM	Giulia Micheli
17/04/2021, 19:19:56	03806999	Patient 03806999 needs psychological support	PSYCHOLOGICAL SUPPORT	Giulia Micheli
17/04/2021, 16:30:02	03806999	Subject 03806999 has «Questionario sulla Salute - EQ-5D-3L» as pending		Subject answered the questionnaire.
15/04/2021, 16:30:01	15929649	Subject 15929649 has «Questionario sulla Salute - EQ-5D-3L» as pending		Giulia Micheli
10/04/2021, 09:10:41	11294295	Patient 11294295 needs psychological support	PSYCHOLOGICAL SUPPORT	Cristina Seguiti

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**DYNAMIC LISTS**

Drag a column header here to group by that column

Search...

Title	Description	Created On	Status	Actions
test		05/03/2021	ACTIVE	

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Patient Criteria

Status: Enrolled Gender: ---Select--- Age:

Baseline Date: Termination Date: Tags: Not Synced

Show: 10

Sr.No	Subject ID	Registration Date
No matching records found		

Showing 0 to 0 of 0 entries (filtered from 95 total entries)

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Description	Type	Content	Recipients	Sent On (User Time)	Sent On (Your Time)
<i>Instant Notification</i>	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	<i>Based on each user's time...</i>	02/02/2022, 16:24:58
<i>Instant Notification</i>	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	<i>Based on each user's time...</i>	02/02/2022, 14:21:49
<i>Instant Notification</i>	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	<i>Based on each user's time...</i>	31/01/2022, 16:28:38
<i>Instant Notification</i>	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	<i>Based on each user's time...</i>	31/01/2022, 16:27:03
<i>Instant Notification</i>	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	<i>Based on each user's time...</i>	28/01/2022, 20:46:50
<i>Instant Notification</i>	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	<i>Based on each user's time...</i>	28/01/2022, 12:33:32
<i>Instant Notification</i>	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	<i>Based on each user's time...</i>	27/01/2022, 15:43:10
<i>Instant Notification</i>	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	<i>Based on each user's time...</i>	27/01/2022, 14:39:19
<i>Instant Notification</i>	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	<i>Based on each user's time...</i>	21/01/2022, 15:03:14
<i>Instant Notification</i>	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	<i>Based on each user's time...</i>	18/01/2022, 11:54:26

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**SCHEDULED NOTIFICATION**

Drag a column header here to group by that column

Description	Type	Content	Status	Actions
Weekly COVID19 Self-Assessment	Weekly	<a href="#">View Content</a>	INACTIVE	<a href="#">Edit</a> <a href="#">Delete</a>
Test	Monthly	<a href="#">View Content</a>	INACTIVE	<a href="#">Edit</a> <a href="#">Delete</a>

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**EDIT SCHEDULED NOTIFICATION**

Description: \*  Status:

Schedule Type:  Regular  By Baseline Date

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
Title	Codename	Type	No. of Questions	Types	Created On	Status	Actions
Initial Questionnaire	ISPRINT_COVID19	Simple	11	INITIAL	16/03/2020	ACTIVE	[Icons]
Body temperature	ISPRINT_FEVER	Simple	1	ADD-EVENT-EDIT	18/03/2020	ACTIVE	[Icons]
Cough	ISPRINT_COUGH	Simple	2	ADD-EVENT-EDIT	18/03/2020	ACTIVE	[Icons]
Fatigue	ISPRINT_FATIGUE	Simple	1	ADD-EVENT-EDIT	18/03/2020	ACTIVE	[Icons]
Difficulty of breath...	ISPRINT_DIBREATH	Simple	1	ADD-EVENT-EDIT	26/03/2020	ACTIVE	[Icons]
Headache	ISPRINT_HEADACHE	Simple	1	ADD-EVENT-EDIT	26/03/2020	ACTIVE	[Icons]
Muscle Pain	ISPRINT_MUSCLES	Simple	2	ADD-EVENT-EDIT	26/03/2020	ACTIVE	[Icons]
Abdominal pain	ISPRINT_ABDPAIN	Simple	1	ADD-EVENT-EDIT	26/03/2020	ACTIVE	[Icons]
Diarrhea	ISPRINT_DIARRHEA	Simple	1	ADD-EVENT-EDIT	26/03/2020	ACTIVE	[Icons]
Nausea	ISPRINT_NAUSEA	Simple	1	ADD-EVENT-EDIT	26/03/2020	ACTIVE	[Icons]
Loss of Smell/Taste	ISPRINT_TASTESM	Simple	1	ADD-EVENT-EDIT	31/03/2020	ACTIVE	[Icons]
Oxygen saturation	ISPRINT_O2KIGEN	Simple	1	ADD-EVENT-EDIT	07/04/2020	ACTIVE	[Icons]
Weekly Questionnaire...	ISPRINT_COVID19_D	Simple	4		08/04/2020	ACTIVE	[Icons]
Add or take photo	PHOTO	Simple	1	ADD-EVENT-EDIT	24/04/2020	ACTIVE	[Icons]
Emotional State	EMOSTATUS	Simple	1	ADD-EVENT-EDIT	06/05/2020	INACTIVE	[Icons]

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
Search...

Name	Actions
App Manager	<a href="#">Edit</a> <a href="#">Delete</a>
Application	<a href="#">Edit</a> <a href="#">Delete</a>
CRO Assistant	<a href="#">Edit</a> <a href="#">Delete</a>
Investigator	<a href="#">Edit</a> <a href="#">Delete</a>
Organization Admin	<a href="#">Edit</a> <a href="#">Delete</a>
System Administrator	<a href="#">Edit</a> <a href="#">Delete</a>
test	<a href="#">Edit</a> <a href="#">Delete</a>
Viewer	<a href="#">Edit</a> <a href="#">Delete</a>

10 20 50 100 All

**CREATE ROLE**

Name: \*

The Role

Role Access Permissions:

- All
  - Dashboard
    - Export
  - Subject
    - Create
    - Edit
    - Details
    - Personal Details
    - Delete
    - View Email
  - Reported Events
    - Export
    - Details
  - Reported Questionnaires
    - Export
  - Notifications
    - Measured Activity Analytics
  - Consent Actions
    - Export
  - Instant Notification
  - Scheduled Notification
    - Create
    - Edit
    - Delete
    - Details
    - ViewContent
    - Export
  - Questionnaire
    - Create
    - Edit

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The screenshot displays the Healthentia User Management interface. The left sidebar contains navigation options: Symptom Tracker, DASHBOARD, STUDY, SUBJECTS, USERS, User List, User Invitations, ALERTS, ECONSENT, NOTIFICATIONS, QUESTIONNAIRES, SYSTEM, DEVICES, and AUDIT TRAIL. The main content area is titled 'USER MANAGEMENT' and features a table of users with columns for Email Address, Role, Organization, Created On, and Actions. Below the table is a pagination control showing 5 items per page and 2 pages total.

Email Address	Role	Organization	Created On	Actions
8afe7900-8d69-4030-b546-8c94f2165170@healthentia.com	Application	ORG_1	02/10/2020	[Edit] [More]
bb426815-dbe6-4182-862e-56139ef4d9cd@healthentia.com	System Administrator	ORG_1	01/03/2021	[Edit] [More]
ef9be18f-06ec-4bc1-b551-a8fada21ffe8@healthentia.com	System Administrator	ORG_1	02/10/2020	[Edit] [More]
f0c0b242-19b7-49b3-9b56-cb442ae32b95@healthentia.com	System Administrator	ORG_1	02/10/2020	[Edit] [More]
f6340709-2eb5-48d7-b755-38c5db756c91@healthentia.com	System Administrator	ORG_1	02/10/2020	[Edit] [More]

The 'CREATE USER' form includes the following fields:

- First Name: \*
- Last Name: \*
- Email Address: \*
- Password: \*
- Confirm Password: \*
- Role: \* (Dropdown menu)
- Time Zone: \* (Dropdown menu)
- Sites: \*

Buttons for 'Cancel' and 'Submit' are located at the bottom right of the form.

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Email	Action	Component	Message	Time Stamp
e9be18f-06ec-4bc1-b551-a8fada21...	View	Subject	Subject List on Study: Symptom Trac...	2022/02/03 09:26:32
e9be18f-06ec-4bc1-b551-a8fada21...	Login	User	SUCCESS	2022/02/03 09:26:26
6cabf309-039d-4ade-a5c5-875a763...	View	Medication	Medication List of Subject: 5e531aa...	2022/02/03 09:25:39
6cabf309-039d-4ade-a5c5-875a763...	View	PatientQuestionnaire	Subject Questionnaire List of Subjec...	2022/02/03 09:25:39
6cabf309-039d-4ade-a5c5-875a763...	View	EConsent	EConsent List of Subject: 5e531aae-...	2022/02/03 09:25:39
6cabf309-039d-4ade-a5c5-875a763...	View	Alert	Alert List of Subject: 5e531aae-76e4...	2022/02/03 09:25:38
6cabf309-039d-4ade-a5c5-875a763...	View	Conversation	Alert List of Subject: 5e531aae-76e4...	2022/02/03 09:25:38
6cabf309-039d-4ade-a5c5-875a763...	View	PatientQuestionnaire	Subject Questionnaire List of Subjec...	2022/02/03 09:25:38
6cabf309-039d-4ade-a5c5-875a763...	View	Subject	Subject Item on Study: Symptom Tra...	2022/02/03 09:25:35
6cabf309-039d-4ade-a5c5-875a763...	View	Subject	Subject List on Study: Symptom Trac...	2022/02/03 09:25:35

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