NAVODILA ZA UPORABO

ZDRAVSTVENI STROKOVNJAKI

 IFU Version v1

INNOVATION SPRINT SRL
Clos Chapelle-aux-Champs 30, bte 1.30.30
1200 Brussels, Belgium

Document History

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Overall
3 studies
422 participants

Tracked Activity
60739 days
128 subjects

Adherence
2804 events from 239 participants
385/624 answered questionnaires

Last week's activity
1 subjects with reported events
1 participants

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Search... [Icons]

Icon	Name	Organization	Codename	Disease	Patients	Status	Alerts	Actions
	Test Study-I	ORG_1	testirida		0	ONGOING	0	
	Interface	ORG_6	Interface	HIV Remote Patient ...	66	ONGOING	10	
	Symptom Tracker	ORG_1	COVID19		356	ONGOING	285	

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Account Settings

Name: * Therapeutic Area: Codename: *

Protocol: Duration In Months: Start Date: 2021/05/14 End Date:

Primary Language: * Add More Languages: Study Logo:

Choose File No file chosen

Additional Services

Video-Chat Multicenter Interactive Map

MOBILE CONFIGURATION

Study Information Page:

Leave empty for no Information Page Study Logo on App:

Choose File No file chosen

Pop-up for Study's terms & Privacy

Short Term's Description: Leave empty for no Short Terms Description

Terms of Use URL:

Leave empty for no Terms of Use

Privacy Policy URL:

Leave empty for no Privacy Policy

Available Consents

General Terms Allowing Reporting Contacted by doctor in case of need

Additional App Widgets

Treatment Reminder Two-Factor Authentication Liquid consumption Sleep Widget

Activity Widget (Available Integrations)

Fitbit Sync Garmin Sync Apple Health Sync Android Sync

SUBJECT DETAIL PAGE CONFIGURATION

First/Last Name Fiscal Number Study Dates Severity

Managed Fitbit Tracker Description Mobile Description Therapeutic Area

Subject Id Email Phone Number Status

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



The screenshot shows the 'TAG MANAGEMENT' interface. On the left is a navigation menu with 'Tag List' selected. The main area displays a table of tags with the following data:

Title	Created On	Type	Status	Actions
Low	23/03/2020	Auto	ACTIVE	[Edit] [Delete]
Medium	23/03/2020	Auto	ACTIVE	[Edit] [Delete]
High	23/03/2020	Auto	ACTIVE	[Edit] [Delete]
Dev	01/04/2020	Custom	ACTIVE	[Edit] [Delete]
Not Synced	29/05/2020	Auto	ACTIVE	[Edit] [Delete]

At the bottom of the table, there are pagination controls showing '5' items per page, '10' total items, and 'All' filters. A navigation bar at the bottom right shows page numbers 1, 2, 3.

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The screenshot shows the 'CREATE DYNAMIC ALERT' form. The form includes the following fields:

- Name:** A text input field.
- Status:** A dropdown menu currently set to 'Active'.
- Portal Alert Message:** A text input field with a note: '(You can use {shortid} to include recipient's Short Id in your message)'. Below it is a small preview box.
- Push Notification Message:** A text input field.
- Email List:** A text input field with a note: '(Enter emails separated by spaces or commas)'.
- Email Message:** A larger text input field with a note: '(You can use {shortid} to include recipient's Short Id in your message)'. Below it is a small preview box.
- API Endpoint to Call:** A text input field.
- API Alert Codename:** A text input field.
- Tags to Apply:** A text input field.

At the bottom right, there are two buttons: 'Cancel' and 'Create'.

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ALERTS

Filter by Tags: COVID19 POSITIVE, DOSE 1 COVIDVAX, DOSE 2 COVIDVAX, FEEDBACK, FITBIT TOKEN EXPIRED, HIGH COVID RISK, HIVSRQ SYMPTOM, LOW ADHERENCE, LOW COVID RISK, NOT SYNCED, PAZIENTE FRAGILE, PSYCHIATRIC SUPPORT, PSYCHOLOGICAL SUPPORT

Show: 50

Triggered On (Browser Time)	Subject Id	Alert Message	Tag	Manager
21/04/2021, 14:21:27	20575866	Patient 20575866 needs psychological support	PSYCHOLOGICAL SUPPORT	Giulia Micheli
20/04/2021, 16:30:02	14455901	Subject 14455901 has «Questionario sulla Salute - EQ-5D-3L» as pending		Giulia Micheli
18/04/2021, 16:30:01	04588190	Subject 04588190 has «Questionario sulla Salute - EQ-5D-3L» as pending		Giulia Micheli
18/04/2021, 13:15:19	14455901	Patient 14455901 is identified as FRAIL	PAZIENTE FRAGILE	Giulia Micheli
18/04/2021, 13:12:09	14455901	Subject 14455901 has an HIVSRQ symptom alert	HIVSRQ SYMPTOM	Giulia Micheli
17/04/2021, 19:19:56	03806999	Patient 03806999 needs psychological support	PSYCHOLOGICAL SUPPORT	Giulia Micheli
17/04/2021, 16:30:02	03806999	Subject 03806999 has «Questionario sulla Salute - EQ-5D-3L» as pending		Subject answered the questionnaire.
15/04/2021, 16:30:01	15929649	Subject 15929649 has «Questionario sulla Salute - EQ-5D-3L» as pending		Giulia Micheli
10/04/2021, 09:10:41	11294295	Patient 11294295 needs psychological support	PSYCHOLOGICAL SUPPORT	Cristina Seguiti

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DYNAMIC LISTS

Drag a column header here to group by that column

Search...

Title	Description	Created On	Status	Actions
test		05/03/2021	ACTIVE	

5 10 All

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Healthentia
Real World Outcomes

Symptom Tracker

- DASHBOARD
- STUDY
- SUBJECTS
- USERS
- ALERTS
- ECONSENT
- NOTIFICATIONS
- QUESTIONNAIRES
- SYSTEM
- DEVICES
- AUDIT TRAIL

DETAILS DYNAMIC LISTS

Title: Status:

Description:

Patient Criteria

Status: Gender: Age:

Baseline Date: Termination Date: Tags:

Back

Show: 10

Sr.No	Subject ID	Registration Date
No matching records found		

Showing 0 to 0 of 0 entries (filtered from 95 total entries)

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Description	Type	Content	Recipients	Sent On (User Time)	Sent On (Your Time)
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	02/02/2022, 16:24:58
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	02/02/2022, 14:21:49
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	31/01/2022, 16:28:38
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	31/01/2022, 16:27:03
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	28/01/2022, 20:46:50
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	28/01/2022, 12:33:32
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	27/01/2022, 15:43:10
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	27/01/2022, 14:39:19
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	21/01/2022, 15:03:14
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	18/01/2022, 11:54:26

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SCHEDULED NOTIFICATION

Drag a column header here to group by that column

Description	Type	Content	Status	Actions
Weekly COVID19 Self-Assessment	Weekly	View Content	INACTIVE	Edit Delete
Test	Monthly	View Content	INACTIVE	Edit Delete

10 50 All

Slika 20: Načrtovana obvestila

EDIT SCHEDULED NOTIFICATION

Description: * Status:

Schedule Type: Regular By Baseline Date

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
Subject	Title	Last Message	Sent On	Category	Reply
SBCYT	Information	Subject: Like arm pain	04/05/2021, 18:20:28	HEALTH	

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
Title	Version	Files	Questions	Actions
eConsent	1	Video		

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

Title	Codename	Type	No. of Questions	Types	Created On	Status	Actions
Initial Questionnaire	ISPRINT_COVID19	Simple	11	INITIAL	16/03/2020	ACTIVE	[edit] [delete]
Body temperature	ISPRINT_FEVER	Simple	1	ADD-EVENT-LIST	18/03/2020	ACTIVE	[edit] [delete]
Cough	ISPRINT_COUGH	Simple	2	ADD-EVENT-LIST	18/03/2020	ACTIVE	[edit] [delete]
Fatigue	ISPRINT_FATIGUE	Simple	1	ADD-EVENT-LIST	18/03/2020	ACTIVE	[edit] [delete]
Difficulty of breath...	ISPRINT_DBREATH	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	[edit] [delete]
Headache	ISPRINT_HEADACHE	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	[edit] [delete]
Muscle Pain	ISPRINT_MUSCLES	Simple	2	ADD-EVENT-LIST	26/03/2020	ACTIVE	[edit] [delete]
Abdominal pain	ISPRINT_ABDPAIN	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	[edit] [delete]
Diarrhea	ISPRINT_DIARRHEA	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	[edit] [delete]
Nausea	ISPRINT_NAUSEA	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	[edit] [delete]
Loss of Smell/Taste	ISPRINT_TASTESM	Simple	1	ADD-EVENT-LIST	31/03/2020	ACTIVE	[edit] [delete]
Oxygen saturation	ISPRINT_OXYGEN	Simple	1	ADD-EVENT-LIST	07/04/2020	ACTIVE	[edit] [delete]
Weekly Questionna...	ISPRINT_COVID19_D	Simple	4		08/04/2020	ACTIVE	[edit] [delete]
Add or take photo	PHOTO	Simple	1	ADD-EVENT-LIST	24/04/2020	ACTIVE	[edit] [delete]
Emotional State	EMOSTATUS	Simple	1	ADD-EVENT-LIST	06/05/2020	INACTIVE	[edit] [delete]

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ROLES & PERMISSIONS FOR ORG...1

Drag a column header here to group by that column

Q Search... [Icons]

Name	Search	Actions
App Manager		[Edit] [Delete]
Application		[Edit] [Delete]
CRO Assistant		[Edit] [Delete]
Investigator		[Edit] [Delete]
Organization Admin		[Edit] [Delete]
System Administrator		[Edit] [Delete]
test		[Edit] [Delete]
Viewer		[Edit] [Delete]

10 20 50 100 All

CREATE ROLE

Name: *

The Role

Role Access Permissions:

- All
 - Dashboard
 - Export
 - Subject
 - Create
 - Edit
 - Details
 - Personal Details
 - Delete
 - View Email
 - Reported Events
 - Export
 - Details
 - Reported Questionnaires
 - Export
 - Notifications
 - Measured Activity Analytics
 - Consent Actions
 - Export
 - Instant Notification
 - Scheduled Notification
 - Create
 - Edit
 - Delete
 - Details
 - ViewContent
 - Export
 - Questionnaire
 - Create
 - Edit

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The screenshot displays the Healthentia web portal interface, specifically the 'USER MANAGEMENT' section. The left sidebar contains navigation options: Symptom Tracker, DASHBOARD, STUDY, SUBJECTS, USERS, User List, User Invitations, ALERTS, ECONSENT, NOTIFICATIONS, QUESTIONNAIRES, SYSTEM, DEVICES, and AUDIT TRAIL. The main content area shows a table of users with columns for Email Address, Role, Organization, Created On, and Actions. Below the table is a pagination control showing 5 items per page and 2 pages total.

Email Address	Role	Organization	Created On	Actions
8afe7900-8d69-4030-b546-8c94f2165170@healthentia.com	Application	ORG_1	02/10/2020	[Edit] [More]
bb426815-dbe6-4182-862e-56139ef4d9cd@healthentia.com	System Administrator	ORG_1	01/03/2021	[Edit] [More]
ef9be18f-06ec-4bc1-b551-a8fada21ffe8@healthentia.com	System Administrator	ORG_1	02/10/2020	[Edit] [More]
f0c0b242-19b7-49b3-9b56-cb442ae32b95@healthentia.com	System Administrator	ORG_1	02/10/2020	[Edit] [More]
f8340709-2eb5-48d7-b755-38c5db756c91@healthentia.com	System Administrator	ORG_1	02/10/2020	[Edit] [More]

Below the table, there is a 'CREATE USER' form with the following fields:

- First Name: *
- Last Name: *
- Email Address: *
- Password: *
- Confirm Password: *
- Role: *
- Time Zone: *
- Sites:

At the bottom right of the form, there are 'Cancel' and 'Submit' buttons.

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PORTAL LOG

Drag a column header here to group by that column

Email	Action	Component	Message	Time Stamp
e9be18f-06ec-4bc1-b551-a8fada21...	View	Subject	Subject List on Study: Symptom Trac...	2022/02/03 09:26:32
e9be18f-06ec-4bc1-b551-a8fada21...	Login	User	SUCCESS	2022/02/03 09:26:26
6cabf309-039d-4ade-a5c5-875a763...	View	Medication	Medication List of Subject: 5e531aa...	2022/02/03 09:25:39
6cabf309-039d-4ade-a5c5-875a763...	View	PatientQuestionnaire	Subject Questionnaire List of Subjec...	2022/02/03 09:25:39
6cabf309-039d-4ade-a5c5-875a763...	View	EConsent	EConsent List of Subject: 5e531aae-...	2022/02/03 09:25:39
6cabf309-039d-4ade-a5c5-875a763...	View	Alert	Alert List of Subject: 5e531aae-76e4...	2022/02/03 09:25:38
6cabf309-039d-4ade-a5c5-875a763...	View	Conversation	Alert List of Subject: 5e531aae-76e4...	2022/02/03 09:25:38
6cabf309-039d-4ade-a5c5-875a763...	View	PatientQuestionnaire	Subject Questionnaire List of Subjec...	2022/02/03 09:25:38
6cabf309-039d-4ade-a5c5-875a763...	View	Subject	Subject Item on Study: Symptom Tra...	2022/02/03 09:25:35
6cabf309-039d-4ade-a5c5-875a763...	View	Subject	Subject List on Study: Symptom Trac...	2022/02/03 09:25:35

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