NAUDOJIMO INSTRUKCIJA

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Icon	Name	Organization	Codename	Disease	Patients	Status	Alerts	Actions
	Test Study-I	ORG_1	testirida		0	ONGOING	0	
	Interface	ORG_6	Interface	HIV Remote Patient ...	66	ONGOING	10	
	Symptom Tracker	ORG_1	COVID19		356	ONGOING	285	

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Account Settings

Name: * Therapeutic Area: Codename: *

Protocol: Duration In Months: Start Date: 2021/05/14 End Date:

Primary Language: * Add More Languages: Study Logo:

Choose File No file chosen

Additional Services

Video-Chat Multicenter Interactive Map

MOBILE CONFIGURATION

Study Information Page: Leave empty for no Information Page Study Logo on App:

Choose File No file chosen

Pop-up for Study's terms & Privacy

Short Term's Description: Leave empty for no Short Terms Description Terms of Use URL:

Leave empty for no Terms of Use

Privacy Policy URL: Leave empty for no Privacy Policy

Available Consents

General Terms Allowing Reporting Contacted by doctor in case of need

Additional App Widgets

Treatment Reminder Two-Factor Authentication Liquid consumption Sleep Widget

Activity Widget (Available Integrations)

Fitbit Sync Garmin Sync Apple Health Sync Android Sync

SUBJECT DETAIL PAGE CONFIGURATION

First/Last Name Fiscal Number Study Dates Severity

Managed Fitbit Tracker Description Mobile Description Therapeutic Area

Subject Id Email Phone Number Status

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3.1 INTERACTIVE MAP

At the configuration of the study, the Study Admin can select to enable the interactive map if there are questions regarding location. Portal users can select parameters through several filters and have an overview of subject outcomes at a geographical level, as presented in Figure 6.




Figure 6: Interactive Map

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

Similarly, to the Email invitation above, the investigator can invite patients by providing them a code. Subjects can enter this code during registration. This is similar to an email invitation, but instead of an email, the investigator can invite patients by providing them with a code. Subjects can enter this code during registration, even if they did not receive an email. Registered subjects are automatically assigned to the specific study. Codes can be generated in bulk – one code for each subject or one common code for all study subjects (see 10 page).

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

The screenshot shows the 'TAG MANAGEMENT' section of the Healthentia portal. The left sidebar contains navigation options: Symptom Tracker, DASHBOARD, STUDY, SUBJECTS, Subject List, Tag List (selected), Dynamic Lists, Custom Fields, Subject Invitations, and Codes. The main area displays a table of tags with the following data:

Title	Created On	Type	Status	Actions
Low	23/03/2020	Auto	ACTIVE	[Edit] [Delete]
Medium	23/03/2020	Auto	ACTIVE	[Edit] [Delete]
High	23/03/2020	Auto	ACTIVE	[Edit] [Delete]
Dev	01/04/2020	Custom	ACTIVE	[Edit] [Delete]
Not Synced	29/05/2020	Auto	ACTIVE	[Edit] [Delete]

At the bottom of the table, there are pagination controls showing '5' items per page and 'All' items. The top right of the interface shows a user profile for 'Innovation Sprint'.

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The screenshot shows the 'CREATE DYNAMIC ALERT' form in the Healthentia portal. The left sidebar is similar to the previous screenshot, with 'Alerts' selected under the 'SUBJECTS' category. The form fields are as follows:

- Name:** * (text input)
- Status:** (dropdown menu, currently set to 'Active')
- Portal Alert Message:** (text input, with a note: "(You can use [shortid] to include recipient's Short Id in your message)")
- Push Notification Message:** (text input)
- Email List:** (text input, with a note: "(Enter emails separated by spaces or commas)")
- Email Message:** (text input, with a note: "(You can use [shortid] to include recipient's Short Id in your message)")
- API Endpoint to Call:** (text input)
- API Alert Codename:** (text input)
- Tags to Apply:** (text input)

At the bottom right of the form, there are 'Cancel' and 'Create' buttons.

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Triggered On (Browser Time)	Subject Id	Alert Message	Tag	Manager
21/04/2021, 14:21:27	20575866	Patient 20575866 needs psychological support	PSYCHOLOGICAL SUPPORT	Giulia Michelli
20/04/2021, 16:30:02	14455901	Subject 14455901 has «Questionario sulla Salute - EQ-5D-3L» as pending		Giulia Michelli
18/04/2021, 16:30:01	04588190	Subject 04588190 has «Questionario sulla Salute - EQ-5D-3L» as pending		Giulia Michelli
18/04/2021, 13:15:19	14455901	Patient 14455901 is identified as FRAIL	PAZIENTE FRAGILE	Giulia Michelli
18/04/2021, 13:12:09	14455901	Subject 14455901 has an HIVSRQ symptom alert	HIVSRQ SYMPTOM	Giulia Michelli
17/04/2021, 19:19:56	03806999	Patient 03806999 needs psychological support	PSYCHOLOGICAL SUPPORT	Giulia Michelli
17/04/2021, 16:30:02	03806999	Subject 03806999 has «Questionario sulla Salute - EQ-5D-3L» as pending		Subject answered the questionnaire.
15/04/2021, 16:30:01	15929649	Subject 15929649 has «Questionario sulla Salute - EQ-5D-3L» as pending		Giulia Michelli
10/04/2021, 09:10:41	11294295	Patient 11294295 needs psychological support	PSYCHOLOGICAL SUPPORT	Cristina Seguiti

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Title	Description	Created On	Status	Actions
test	test	05/03/2021	ACTIVE	[Edit] [Delete]

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Healthentia
Real World Outcomes

Symptom Tracker

- DASHBOARD
- STUDY
- SUBJECTS
- USERS
- ALERTS
- ECONSENT
- NOTIFICATIONS
- QUESTIONNAIRES
- SYSTEM
- DEVICES
- AUDIT TRAIL

DETAILS DYNAMIC LISTS

Title: Status:

Description:

Patient Criteria

Status: Gender: Age:

Baseline Date: Termination Date: Tags:

Back

Show: 10

Sr.No	Subject ID	Registration Date
No matching records found		

Showing 0 to 0 of 0 entries (filtered from 95 total entries)

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Description	Type	Content	Recipients	Sent On (User Time)	Sent On (Your Time)
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	02/02/2022, 16:24:58
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	02/02/2022, 14:21:49
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	31/01/2022, 16:28:38
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	31/01/2022, 16:27:03
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	28/01/2022, 20:46:50
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	28/01/2022, 12:33:32
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	27/01/2022, 15:43:10
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	27/01/2022, 14:39:19
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	21/01/2022, 15:03:14
< >Instant Notification</ >	Instant	View Content	View Recipients (1)	< >Based on each user's time...	18/01/2022, 11:54:26

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SCHEDULED NOTIFICATION

Drag a column header here to group by that column

Description	Type	Content	Status	Actions
Weekly COVID19 Self-Assessment	Weekly	View Content	INACTIVE	Edit Delete
Test	Monthly	View Content	INACTIVE	Edit Delete

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Description: *
 Test

Status:
 Inactive

Schedule Type:
 Regular By Baseline Date

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Subject	Title	Last Message	Sent On	Category	Reply
SBCYT	Information	Subject: Like arm pain	04/05/2021, 18:20:28	HEALTH	

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CONVERSATION WITH KONSTANTINA KOSTOPOULOU: ΔΩΝΒΓ

START TELECONSULTATION

Apr 22, 2021

Γιωργι
Konstantina Kostopoulou 12:53 PM

hello
12:54 PM

Teleconsultation
Teleconsultation request was sent.
12:55 PM

Teleconsultation has been completed.
12:57 PM

Teleconsultation
Teleconsultation request was sent.
12:58 PM

Teleconsultation has been completed.
12:58 PM

Teleconsultation
Teleconsultation request was sent.
6:57 PM

Teleconsultation has been completed.
6:58 PM

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

Title	Codename	Type	No. of Questions	Types	Created On	Status	Actions
Initial Questionnaire	ISPRINT_COVID19	Simple	11	INITIAL	16/03/2020	ACTIVE	
Body temperature	ISPRINT_FEVER	Simple	1	ADD-EVENT-LIST	18/03/2020	ACTIVE	
Cough	ISPRINT_COUGH	Simple	2	ADD-EVENT-LIST	18/03/2020	ACTIVE	
Fatigue	ISPRINT_FATIGUE	Simple	1	ADD-EVENT-LIST	18/03/2020	ACTIVE	
Difficulty of breath...	ISPRINT_DBREATH	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	
Headache	ISPRINT_HEADACHE	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	
Muscle Pain	ISPRINT_MUSCLES	Simple	2	ADD-EVENT-LIST	26/03/2020	ACTIVE	
Abdominal pain	ISPRINT_ABDPAIN	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	
Diarrhea	ISPRINT_DIARRHEA	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	
Nausea	ISPRINT_NAUSEA	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	
Loss of Smell/Taste	ISPRINT_TASTESM	Simple	1	ADD-EVENT-LIST	31/03/2020	ACTIVE	
Oxygen saturation	ISPRINT_OXYGEN	Simple	1	ADD-EVENT-LIST	07/04/2020	ACTIVE	
Weekly Questionna...	ISPRINT_COVID19_D	Simple	4		08/04/2020	ACTIVE	
Add or take photo	PHOTO	Simple	1	ADD-EVENT-LIST	24/04/2020	ACTIVE	
Emotional State	EMOSTATUS	Simple	1	ADD-EVENT-LIST	06/05/2020	INACTIVE	

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The screenshot shows the 'VIEW QUESTIONNAIRE' page in the Healthentia system. The left sidebar contains navigation options: Symptom Tracker, DASHBOARD, STUDY, SUBJECTS, USERS, ALERTS, COMMENT, NOTIFICATIONS, QUESTIONNAIRES, SYSTEM, DEVICES, and AUDIT TRAIL. The main content area is titled 'VIEW QUESTIONNAIRE' and displays details for a questionnaire with the title 'Body temperature' and codename 'ISPRINT_FEVER'. The status is 'Active'. The 'Types' field is set to 'Add-Event-List'. There is an 'Image' field with a placeholder, a 'Color' field set to 'E58853' (a red color swatch), and a 'Description' field. Below this, a 'QUESTION LIST' section shows a table with one entry:

No	Code	Question	UI Control	Created On
1	FEVER	What is your body temperature?	Numerical Input	18/03/2020

At the bottom, it indicates 'Showing 1 to 1 of 1 entries' and has a 'Cancel' button.

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The screenshot shows the 'VIEW QUESTIONNAIRE' page for an 'Initial Questionnaire' with codename 'ISPRINT_COVID19' and status 'Active'. The 'Types' field is set to 'Initial'. There is an 'Image' field with a placeholder, a 'Color' field set to '808080' (a dark grey color swatch), and a 'Description' field. Below this, a 'QUESTION LIST' section shows a table with five entries:

No	Code	Question	UI Control	Created On
1	BIRTHDT	Enter your Birthdate	Date	13/04/2020
2	GENDER	Define your gender:	Single Choice	16/03/2020
3	DIAGNOSE	Have you been diagnosed with COVID-19?	Single Choice	16/03/2020
4	SYMPTOMS	Do you currently have any of the following symptoms?	Multiple Choice	16/03/2020
5	MEDICIN	Do you usually get any of these medications?	Single Choice	01/04/2020

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Q Search... [Icons]

Name	Actions
App Manager	[Edit] [Delete]
Application	[Edit] [Delete]
CRO Assistant	[Edit] [Delete]
Investigator	[Edit] [Delete]
Organization Admin	[Edit] [Delete]
System Administrator	[Edit] [Delete]
test	[Edit] [Delete]
Viewer	[Edit] [Delete]

10 20 50 100 All

CREATE ROLE

Name: *
The Role


Role Access Permissions:

- All
 - Dashboard
 - Export
 - Subject
 - Create
 - Edit
 - Details
 - Personal Details
 - Delete
 - View Email
 - Reported Events
 - Export
 - Details
 - Reported Questionnaires
 - Export
 - Notifications
 - Measured Activity Analytics
 - Consent Actions
 - Export
 - Instant Notification
 - Scheduled Notification
 - Create
 - Edit
 - Delete
 - Details
 - ViewContent
 - Export
 - Questionnaire
 - Create
 - Edit

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
The screenshot shows the 'USER MANAGEMENT' section of the Healthentia web portal. The interface includes a sidebar with navigation options like 'Symptom Tracker', 'DASHBOARD', 'STUDY', 'SUBJECTS', 'USERS', and 'User List'. The main content area displays a table of users with columns for 'Email Address', 'Role', 'Organization', 'Created On', and 'Actions'. There are search bars for each column and a 'Search...' field. The table lists five users, all with the role 'System Administrator' and organization 'ORG_1'. The 'Created On' dates range from 02/10/2020 to 01/03/2021. At the bottom, there are pagination controls showing '5' and '10' items per page, and a page number '2'.

Email Address	Role	Organization	Created On	Actions
8afe7900-8d69-4030-b546-8c94f2165170@healthentia.com	Application	ORG_1	02/10/2020	[Edit] [More]
bb426815-dbe6-4182-862e-56139ef4d9cd@healthentia.com	System Administrator	ORG_1	01/03/2021	[Edit] [More]
ef9be18f-06ec-4bc1-b551-a8fada21ffe8@healthentia.com	System Administrator	ORG_1	02/10/2020	[Edit] [More]
f0c0b242-19b7-49b3-9b56-cb442ae32b95@healthentia.com	System Administrator	ORG_1	02/10/2020	[Edit] [More]
f8340709-2eb5-48d7-b755-38c5db756c91@healthentia.com	System Administrator	ORG_1	02/10/2020	[Edit] [More]

The screenshot shows the 'CREATE USER' form in the Healthentia web portal. The form includes input fields for 'First Name', 'Last Name', 'Email Address', 'Password', and 'Confirm Password'. There are also dropdown menus for 'Role' and 'Time Zone', and a text area for 'Sites'. The form is styled with a purple sidebar and a white main content area. At the bottom right, there are 'Cancel' and 'Submit' buttons.

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Email	Action	Component	Message	Time Stamp
ef9be18f-06ec-4bc1-b551-a8fada21...	View	Subject	Subject List on Study: Symptom Trac...	2022/02/03 09:26:32
ef9be18f-06ec-4bc1-b551-a8fada21...	Login	User	SUCCESS	2022/02/03 09:26:26
6cabf309-039d-4ade-a5c5-875a763...	View	Medication	Medication List of Subject: 5e531aa...	2022/02/03 09:25:39
6cabf309-039d-4ade-a5c5-875a763...	View	PatientQuestionnaire	Subject Questionnaire List of Subjec...	2022/02/03 09:25:39
6cabf309-039d-4ade-a5c5-875a763...	View	EConsent	EConsent List of Subject: 5e531aae-...	2022/02/03 09:25:39
6cabf309-039d-4ade-a5c5-875a763...	View	Alert	Alert List of Subject: 5e531aae-76e4...	2022/02/03 09:25:38
6cabf309-039d-4ade-a5c5-875a763...	View	Conversation	Alert List of Subject: 5e531aae-76e4...	2022/02/03 09:25:38
6cabf309-039d-4ade-a5c5-875a763...	View	PatientQuestionnaire	Subject Questionnaire List of Subjec...	2022/02/03 09:25:38
6cabf309-039d-4ade-a5c5-875a763...	View	Subject	Subject Item on Study: Symptom Tra...	2022/02/03 09:25:35
6cabf309-039d-4ade-a5c5-875a763...	View	Subject	Subject List on Study: Symptom Trac...	2022/02/03 09:25:35

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