# NAVODILA ZA UPORABO

## ZDRAVSTVENI STROKOVNJAKI

 IFU Version v2

INNOVATION SPRINT SRL  
Clos Chapelle-aux-Champs 30, bte 1.30.30  
1200 Brussels, Belgium

## Document History

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
Copyright © 2022. All Rights Reserved - [Healthentia](#)

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**Overall**  
3 studies  
422 participants

**Tracked Activity**  
60739 days  
128 subjects

**Adherence**  
2804 events from 239 participants  
385/624 answered questionnaires

**Last week's activity**  
1 subjects with reported events  
1 participants

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Search... [Icons]

Icon	Name	Organization	Codename	Disease	Patients	Status	Alerts	Actions
[Folder Icon]	Test Study-I	ORG_1	testirida		0	ONGOING	0	[Edit] [More]
[Globe Icon]	Interface	ORG_6	Interface	HIV Remote Patient ...	66	ONGOING	10	[Edit] [More]
[Grid Icon]	Symptom Tracker	ORG_1	COVID19		356	ONGOING	285	[Edit] [More]

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##### Account Settings

Name: \* Therapeutic Area: Codename: \*

Protocol: Duration In Months: Start Date: 2021/05/14 End Date:

Primary Language: \* Add More Languages: Study Logo:

Choose File No file chosen

##### Additional Services

Video-Chat  Multicenter  Interactive Map

#### MOBILE CONFIGURATION

##### Study Information Page:

Leave empty for no Information Page Study Logo on App:

Choose File No file chosen

##### Pop-up for Study's terms & Privacy

Short Term's Description: Leave empty for no Short Terms Description

Terms of Use URL:

Leave empty for no Terms of Use

Privacy Policy URL: Leave empty for no Privacy Policy

##### Available Consents

General Terms  Allowing Reporting  Contacted by doctor in case of need

##### Additional App Widgets

Treatment Reminder  Two-Factor Authentication  Liquid consumption  Sleep Widget

##### Activity Widget (Available Integrations)

Fitbit Sync  Garmin Sync  Apple Health Sync  Android Sync

#### SUBJECT DETAIL PAGE CONFIGURATION

First/Last Name  Fiscal Number  Study Dates  Severity

Managed Fitbit  Tracker Description  Mobile Description  Therapeutic Area

Subject Id  Email  Phone Number  Status

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



The screenshot shows the 'TAG MANAGEMENT' interface. On the left is a navigation menu with 'Tag List' selected. The main area displays a table of tags with columns for Title, Created On, Type, Status, and Actions. The table contains five rows of data.

Title	Created On	Type	Status	Actions
Low	23/03/2020	Auto	ACTIVE	[Edit] [Delete]
Medium	23/03/2020	Auto	ACTIVE	[Edit] [Delete]
High	23/03/2020	Auto	ACTIVE	[Edit] [Delete]
Dev	01/04/2020	Custom	ACTIVE	[Edit] [Delete]
Not Synced	29/05/2020	Auto	ACTIVE	[Edit] [Delete]

At the bottom of the table, there are pagination controls showing '5' items per page, '10' total items, and 'All' filter options. A navigation bar at the bottom right shows page numbers 1, 2, 3, and a right arrow.

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The screenshot shows the 'CREATE DYNAMIC ALERT' form. It includes several input fields and a dropdown menu for 'Status'. The form is designed for creating a dynamic alert with specific messages and tags.

Fields include:

- Name: \*
- Status: Active (dropdown)
- Portal Alert Message: (You can use **{shortid}** to include recipient's Short Id in your message)
- Push Notification Message:
- Email List: (Enter emails separated by spaces or commas)
- Email Message: (You can use **{shortid}** to include recipient's Short Id in your message)
- API Endpoint to Call:
- API Alert Codename:
- Tags to Apply:

At the bottom right, there are 'Cancel' and 'Create' buttons.

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Triggered On (Browser Time)	Subject Id	Alert Message	Tag	Manager
21/04/2021, 14:21:27	20575866	Patient 20575866 needs psychological support	PSYCHOLOGICAL SUPPORT	Giulia Micheli
20/04/2021, 16:30:02	14455901	Subject 14455901 has «Questionario sulla Salute - EQ-5D-3L» as pending		Giulia Micheli
18/04/2021, 16:30:01	04588190	Subject 04588190 has «Questionario sulla Salute - EQ-5D-3L» as pending		Giulia Micheli
18/04/2021, 13:15:19	14455901	Patient 14455901 is identified as FRAIL	PAZIENTE FRAGILE	Giulia Micheli
18/04/2021, 13:12:09	14455901	Subject 14455901 has an HIVSRQ symptom alert	HIVSRQ SYMPTOM	Giulia Micheli
17/04/2021, 19:19:56	03806999	Patient 03806999 needs psychological support	PSYCHOLOGICAL SUPPORT	Giulia Micheli
17/04/2021, 16:30:02	03806999	Subject 03806999 has «Questionario sulla Salute - EQ-5D-3L» as pending		Subject answered the questionnaire.
15/04/2021, 16:30:01	15929649	Subject 15929649 has «Questionario sulla Salute - EQ-5D-3L» as pending		Giulia Micheli
10/04/2021, 09:10:41	11294295	Patient 11294295 needs psychological support	PSYCHOLOGICAL SUPPORT	Cristina Seguiti

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Title	Description	Created On	Status	Actions
test	test	05/03/2021	ACTIVE	[Edit] [Delete]

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Healthentia  
Real World Outcomes

Symptom Tracker

DASHBOARD

STUDY

SUBJECTS

USERS

ALERTS

ECONSENT

NOTIFICATIONS

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SYSTEM

DEVICES

AUDIT TRAIL

DETAILS DYNAMIC LISTS

Title:  Status:

Description:

Patient Criteria

Status:  Gender:  Age:

Baseline Date:  Termination Date:  Tags:

Back

Show: 10

Sr.No	Subject ID	Registration Date
No matching records found		

Showing 0 to 0 of 0 entries (filtered from 95 total entries)

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Description	Type	Content	Recipients	Sent On (User Time)	Sent On (Your Time)
< >Instant Notification</ >	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	< >Based on each user's time...	02/02/2022, 16:24:58
< >Instant Notification</ >	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	< >Based on each user's time...	02/02/2022, 14:21:49
< >Instant Notification</ >	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	< >Based on each user's time...	31/01/2022, 16:28:38
< >Instant Notification</ >	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	< >Based on each user's time...	31/01/2022, 16:27:03
< >Instant Notification</ >	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	< >Based on each user's time...	28/01/2022, 20:46:50
< >Instant Notification</ >	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	< >Based on each user's time...	28/01/2022, 12:33:32
< >Instant Notification</ >	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	< >Based on each user's time...	27/01/2022, 15:43:10
< >Instant Notification</ >	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	< >Based on each user's time...	27/01/2022, 14:39:19
< >Instant Notification</ >	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	< >Based on each user's time...	21/01/2022, 15:03:14
< >Instant Notification</ >	Instant	<a href="#">View Content</a>	<a href="#">View Recipients (1)</a>	< >Based on each user's time...	18/01/2022, 11:54:26

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**SCHEDULED NOTIFICATION**

Drag a column header here to group by that column

Description	Type	Content	Status	Actions
Weekly COVID19 Self-Assessment	Weekly	<a href="#">View Content</a>	INACTIVE	<a href="#">Edit</a> <a href="#">Delete</a>
Test	Monthly	<a href="#">View Content</a>	INACTIVE	<a href="#">Edit</a> <a href="#">Delete</a>

10 50 All

Slika 20: Načrtovana obvestila

**EDIT SCHEDULED NOTIFICATION**

Description: \*  
 Test

Status:  
 Inactive

Schedule Type:  
 Regular  By Baseline Date

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
Subject	Title	Last Message	Sent On	Category	Reply
SBCYT	Information	Subject: Like arm pain	04/05/2021, 18:20:28	HEALTH	

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
Title	Version	Files	Questions	Actions
eConsent	1	Video		

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


Title	Codename	Type	No. of Questions	Types	Created On	Status	Actions
Initial Questionnaire	ISPRINT_COVID19	Simple	11	INITIAL	16/03/2020	ACTIVE	[edit] [delete]
Body temperature	ISPRINT_FEVER	Simple	1	ADD-EVENT-LIST	18/03/2020	ACTIVE	[edit] [delete]
Cough	ISPRINT_COUGH	Simple	2	ADD-EVENT-LIST	18/03/2020	ACTIVE	[edit] [delete]
Fatigue	ISPRINT_FATIGUE	Simple	1	ADD-EVENT-LIST	18/03/2020	ACTIVE	[edit] [delete]
Difficulty of breath...	ISPRINT_DBREATH	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	[edit] [delete]
Headache	ISPRINT_HEADACHE	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	[edit] [delete]
Muscle Pain	ISPRINT_MUSCLES	Simple	2	ADD-EVENT-LIST	26/03/2020	ACTIVE	[edit] [delete]
Abdominal pain	ISPRINT_ABDPAIN	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	[edit] [delete]
Diarrhea	ISPRINT_DIARRHEA	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	[edit] [delete]
Nausea	ISPRINT_NAUSEA	Simple	1	ADD-EVENT-LIST	26/03/2020	ACTIVE	[edit] [delete]
Loss of Smell/Taste	ISPRINT_TASTESM	Simple	1	ADD-EVENT-LIST	31/03/2020	ACTIVE	[edit] [delete]
Oxygen saturation	ISPRINT_OXYGEN	Simple	1	ADD-EVENT-LIST	07/04/2020	ACTIVE	[edit] [delete]
Weekly Questionna...	ISPRINT_COVID19_D	Simple	4		08/04/2020	ACTIVE	[edit] [delete]
Add or take photo	PHOTO	Simple	1	ADD-EVENT-LIST	24/04/2020	ACTIVE	[edit] [delete]
Emotional State	EMOSTATUS	Simple	1	ADD-EVENT-LIST	06/05/2020	INACTIVE	[edit] [delete]

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**ROLES & PERMISSIONS FOR ORG\_1**

Drag a column header here to group by that column

Q Search... [Icons]

Name		Actions
App Manager		[Edit] [Delete]
Application		[Edit] [Delete]
CRO Assistant		[Edit] [Delete]
Investigator		[Edit] [Delete]
Organization Admin		[Edit] [Delete]
System Administrator		[Edit] [Delete]
test		[Edit] [Delete]
Viewer		[Edit] [Delete]

10 20 50 100 All

---

**CREATE ROLE**

Name: \*

The Role

Role Access Permissions:

- All
  - Dashboard
    - Export
  - Subject
    - Create
    - Edit
    - Details
    - Personal Details
    - Delete
    - View Email
  - Reported Events
    - Export
    - Details
  - Reported Questionnaires
    - Export
    - Notifications
    - Measured Activity Analytics
  - Consent Actions
    - Export
  - Instant Notification
  - Scheduled Notification
    - Create
    - Edit
    - Delete
    - Details
    - ViewContent
    - Export
  - Questionnaire
    - Create
    - Edit

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The screenshot displays the Healthentia web portal interface. On the left is a navigation sidebar with options like Symptom Tracker, Dashboard, Study, Subjects, Users, Alerts, Econsent, Notifications, Questionnaires, System, Devices, and Audit Trail. The main content area is divided into two sections: 'USER MANAGEMENT' and 'CREATE USER'.

**USER MANAGEMENT**

Drag a column header here to group by that column

Search... [Icons for user, add, edit, delete, refresh]

Email Address	Role	Organization	Created On	Actions
<input type="text" value="8afe7900-8d69-4030-b546-8c94f2165170@healthentia.com"/>	Application	ORG_1	02/10/2020	[Edit] [More]
<input type="text" value="bb426815-dbe6-4182-862e-56139ef4d9cd@healthentia.com"/>	System Administrator	ORG_1	01/03/2021	[Edit] [More]
<input type="text" value="ef9be18f-06ec-4bc1-b551-a8fada21ffe8@healthentia.com"/>	System Administrator	ORG_1	02/10/2020	[Edit] [More]
<input type="text" value="f0c0b242-19b7-49b3-9b56-cb442ae32b95@healthentia.com"/>	System Administrator	ORG_1	02/10/2020	[Edit] [More]
<input type="text" value="f8340709-2eb5-48d7-b755-38c5db756c91@healthentia.com"/>	System Administrator	ORG_1	02/10/2020	[Edit] [More]

5 | 10 | < 1 2 >

**CREATE USER**

First Name: \*  Last Name: \*

Email Address: \*  Password:  Confirm Password:

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**PORTAL LOG**

Drag a column header here to group by that column

Email	Action	Component	Message	Time Stamp
e9be18f-06ec-4bc1-b551-a8fada21...	View	Subject	Subject List on Study: Symptom Trac...	2022/02/03 09:26:32
e9be18f-06ec-4bc1-b551-a8fada21...	Login	User	SUCCESS	2022/02/03 09:26:26
6cabf309-039d-4ade-a5c5-875a763...	View	Medication	Medication List of Subject: 5e531aa...	2022/02/03 09:25:39
6cabf309-039d-4ade-a5c5-875a763...	View	PatientQuestionnaire	Subject Questionnaire List of Subjec...	2022/02/03 09:25:39
6cabf309-039d-4ade-a5c5-875a763...	View	EConsent	EConsent List of Subject: 5e531aae-...	2022/02/03 09:25:39
6cabf309-039d-4ade-a5c5-875a763...	View	Alert	Alert List of Subject: 5e531aae-76e4...	2022/02/03 09:25:38
6cabf309-039d-4ade-a5c5-875a763...	View	Conversation	Alert List of Subject: 5e531aae-76e4...	2022/02/03 09:25:38
6cabf309-039d-4ade-a5c5-875a763...	View	PatientQuestionnaire	Subject Questionnaire List of Subjec...	2022/02/03 09:25:38
6cabf309-039d-4ade-a5c5-875a763...	View	Subject	Subject Item on Study: Symptom Tra...	2022/02/03 09:25:35
6cabf309-039d-4ade-a5c5-875a763...	View	Subject	Subject List on Study: Symptom Trac...	2022/02/03 09:25:35

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